



OPHTHALMED™

LASER PROBES
FOR ENDO-OCULAR
PHOTOCOAGULATION
TREATMENTS

WWW.OPHTHALMED.COM

(877) 4-25 GAUGE (4284)

In business since 2001, Ophthalmed LLC specializes in designing and manufacturing higher quality endo photocoagulation laser probes.

The range of our products is one of the most extensive on the market, backed up by several patent pending applications, and private labeled to major companies in US and worldwide, including well known laser machine manufacturers and other high profile distributors.


We are committed to serving the needs of the ophthalmic industry, and to providing the latest technologies to the surgeons, including custom designs.

Our current product range includes laser and illuminated laser probes in various tip sizes and configurations, including bending (patented). We connect to most laser and illumination machines.

Our probes are manufactured in accordance with a quality system certified to ISO 13485 and cleared/approved by various regulatory agencies, such as US FDA, Health Canada and the European Community.

BENDING LASER PROBES

Protected by US patent 9,370,447
and other patents and patent applications worldwide



Fully slide button forward to straighten the tip prior to penetration or withdrawal; sliding the button all the way forward straightens the tip while sliding the button backward incrementally curves the tip.

REF	TIP	LASER CONNECTOR	
A207000	BENDING, 20G	ALCON / ZEISS / ELLEX	
A807000	BENDING, 23G		
A307000	BENDING, 25G		
AA07000	BENDING, 27G		
A208000	BENDING, 20G	IRIDEX	
A808000	BENDING, 23G		
A308000	BENDING, 25G		
AA08000	BENDING, 27G		
A203000	BENDING, 20G	DORC (non-USA)	
A803000	BENDING, 23G		
A303000	BENDING, 25G		
AA03000	BENDING, 27G		
A209000	BENDING, 20G	QUANTEL SUPRA/VITRA (DIRECT CONNECTION)	
A809000	BENDING, 23G		
A309000	BENDING, 25G		
AA09000	BENDING, 27G		
A20A000	BENDING, 20G	LUMENIS (non-HGM)	
A80A000	BENDING, 23G		
A30A000	BENDING, 25G		
AA0A000	BENDING, 27G		
A202000	BENDING, 20G	NIDEK (THREADED)	
A802000	BENDING, 23G		
A302000	BENDING, 25G		
AA02000	BENDING, 27G		
A20B000	BENDING, 20G	LIGHTMED	
A80B000	BENDING, 23G		
A30B000	BENDING, 25G		
AA0B000	BENDING, 27G		

ALCON / ZEISS / ELLEX LASER CONNECTOR

ALCON ACCURUS AND CONSTELLATION ILLUMINATION CONNECTOR

- A507200 20 g
- B507203 23 g
- B507204 25 g



Protected by US patent 9,370,447 and other patents and patent applications worldwide



19 GAUGE ILLUMINATION CONNECTOR

- A507100 20 g
- B507103 23 g
- B507104 25 g



B&L STELLARIS PC ILLUMINATION CONNECTOR

- A507500 20 g
- B507503 23 g
- B507504 25 g



ALCON OCUTOME ILLUMINATION CONNECTOR

- A507300 20 g
- B507303 23 g
- B507304 25 g



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

- A507400 20 g
- B507403 23 g
- B507404 25 g



IRIDEX LASER CONNECTOR

ALCON ACCURUS AND CONSTELLATION ILLUMINATION CONNECTOR

- A508200 20 g
- B508203 23 g
- B508204 25 g



Protected by US patent 9,370,447 and other patents and patent applications worldwide



19 GAUGE ILLUMINATION CONNECTOR

- A508100 20 g
- B508103 23 g
- B508104 25 g



B&L STELLARIS PC ILLUMINATION CONNECTOR

- A508500 20 g
- B508503 23 g
- B508504 25 g



ALCON OCUTOME ILLUMINATION CONNECTOR

- A508300 20 g
- B508303 23 g
- B508304 25 g



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

- A508400 20 g
- B508403 23 g
- B508404 25 g



DORC (NON-USA) LASER CONNECTOR

ALCON ACCURUS AND CONSTELLATION ILLUMINATION CONNECTOR

A503200 20 g
B503203 23 g
B503204 25 g



Protected by US patent 9,370,447 and other patents and patent applications worldwide



19 GAUGE ILLUMINATION CONNECTOR

A503100 20 g
B503103 23 g
B503104 25 g



B&L STELLARIS PC ILLUMINATION CONNECTOR

A503500 20 g
B503503 23 g
B503504 25 g



ALCON OCUTOME ILLUMINATION CONNECTOR

A503300 20 g
B503303 23 g
B503304 25 g



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

A503400 20 g
B503403 23 g
B503404 25 g



QUANTEL VITRA / SUPRA LASER CONNECTOR (DIRECT)

ALCON ACCURUS AND CONSTELLATION ILLUMINATION CONNECTOR

A509200 20 g
B509203 23 g
B509204 25 g



Protected by US patent 9,370,447 and other patents and patent applications worldwide



19 GAUGE ILLUMINATION CONNECTOR

A509100 20 g
B509103 23 g
B509104 25 g



B&L STELLARIS PC ILLUMINATION CONNECTOR

A509500 20 g
B509503 23 g
B509504 25 g



ALCON OCUTOME ILLUMINATION CONNECTOR

A509300 20 g
B509303 23 g
B509304 25 g



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

A509400 20 g
B509403 23 g
B509404 25 g



LUMENIS (NON-HGM) LASER CONNECTOR

ALCON ACCURUS AND CONSTELLATION ILLUMINATION CONNECTOR

- A50A200 20 g
- B50A203 23 g
- B50A204 25 g



Protected by US patent 9,370,447 and other patents and patent applications worldwide



19 GAUGE ILLUMINATION CONNECTOR

- A50A100 20 g
- B50A103 23 g
- B50A104 25 g



B&L STELLARIS PC ILLUMINATION CONNECTOR

- A50A500 20 g
- B50A503 23 g
- B50A504 25 g



ALCON OCUTOME ILLUMINATION CONNECTOR

- A50A300 20 g
- B50A303 23 g
- B50A304 25 g



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

- A50A400 20 g
- B50A403 23 g
- B50A404 25 g



NIDEK LASER CONNECTOR - THREADED

ALCON ACCURUS AND CONSTELLATION ILLUMINATION CONNECTOR

- A502200 20 g
- B502203 23 g
- B502204 25 g



Protected by US patent 9,370,447 and other patents and patent applications worldwide



19 GAUGE ILLUMINATION CONNECTOR

- A502100 20 g
- B502103 23 g
- B502104 25 g



B&L STELLARIS PC ILLUMINATION CONNECTOR

- A502500 20 g
- B502503 23 g
- B502504 25 g



ALCON OCUTOME ILLUMINATION CONNECTOR

- A502300 20 g
- B502303 23 g
- B502304 25 g



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

- A502400 20 g
- B502403 23 g
- B502404 25 g



LIGHTMED LASER CONNECTOR

ALCON ACCURUS AND CONSTELLATION ILLUMINATION CONNECTOR

- A50B200 20 g
- B50B203 23 g
- B50B204 25 g



Protected by US patent 9,370,447 and other patents and patent applications worldwide



19 GAUGE ILLUMINATION CONNECTOR

- A50B100 20 g
- B50B103 23 g
- B50B104 25 g



B&L STELLARIS PC ILLUMINATION CONNECTOR

- A50B500 20 g
- B50B503 23 g
- B50B504 25 g



ALCON OCUTOME ILLUMINATION CONNECTOR


- A50B300 20 g
- B50B303 23 g
- B50B304 25 g



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

- A50B400 20 g
- B50B403 23 g
- B50B404 25 g





ALCON / ZEISS / ELLEX CONNECTOR

A217000	Straight	
A237000	Angled	20g
A257000	Curved	
A827000	Straight	23g
A877000	Angled	
A327000	Straight	
A377000	Angled	25g
A367000	Curved	
AA27000	Straight	27g
AA67000	Curved	



IRIDEX CONNECTOR

A218000	Straight	
A238000	Angled	20g
A258000	Curved	
A828000	Straight	23g
A878000	Angled	
A328000	Straight	
A378000	Angled	25g
A368000	Curved	
AA28000	Straight	27g
AA68000	Curved	



DORC (NON-US) CONNECTOR

A213000	Straight	
A233000	Angled	20g
A253000	Curved	
A823000	Straight	23g
A873000	Angled	
A323000	Straight	
A373000	Angled	25g
A363000	Curved	
AA23000	Straight	27g
AA63000	Curved	



QUANTEL VITRA / SUPRA CONNECTOR (DIRECT)

A219000	Straight	
A239000	Angled	20g
A259000	Curved	
A829000	Straight	23g
A879000	Angled	
A329000	Straight	
A379000	Angled	25g
A369000	Curved	
AA29000	Straight	27g
AA69000	Curved	



LUMENIS (NON-HGM) CONNECTOR

A21A000	Straight	
A23A000	Angled	20g
A25A000	Curved	
A82A000	Straight	23g
A87A000	Angled	
A32A000	Straight	
A37A000	Angled	25g
A36A000	Curved	
AA2A000	Straight	27g
AA6A000	Curved	



NIDEK CONNECTOR - THREADED


A212000	Straight	
A232000	Angled	20g
A252000	Curved	
A822000	Straight	23g
A872000	Angled	
A322000	Straight	
A372000	Angled	25g
A362000	Curved	
AA22000	Straight	27g
AA62000	Curved	



ADDITIONAL VARIATIONS ARE AVAILABLE.



TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS



NIDEK CONNECTOR - PUSH-IN

A215000	Straight	
A235000	Angled	20g
A255000	Curved	
A825000	Straight	23g
A875000	Angled	
A325000	Straight	
A375000	Angled	25g
A365000	Curved	
AA25000	Straight	27g
AA65000	Curved	



LIGHTMED CONNECTOR

A21b000	Straight	
A23B000	Angled	20g
A25B000	Curved	
A82B000	Straight	23g
A87B000	Angled	
A32B000	Straight	
A37B000	Angled	25g
A36B000	Curved	
AA2B000	Straight	27g
AA6B000	Curved	




HGM CONNECTOR - NO BAR CODE

A214000	Straight	
A234000	Angled	20g
A254000	Curved	
A824000	Straight	23g
A874000	Angled	
A324000	Straight	
A374000	Angled	25g
A364000	Curved	
AA24000	Straight	27g
AA64000	Curved	



STANDARD NON-EXTENDED TIPS



EXTENDED TIPS - TYPICAL



Protected by US patent 9,335,455 and other patents and patent applications worldwide

ALCON / ZEISS / ELLEX LASER CONNECTOR

ALCON ACCURUS AND CONSTELLATION ILLUMINATION CONNECTOR

A527200	Straight	
A547200	Angled	
A567200	Curved	20g
A587200	Straight, Extended	
A597200	Angled, Extended	
A5A7200	Curved, Extended	
B527203	Straight	
B577203	Angled	23g
B587203	Straight, Extended	
B5B7203	Angled, Extended	
B527204	Straight	
B567204	Curved	25g
B587204	Straight, Extended	
B5A7204	Curved, Extended	
B527205	Straight	27g
B567205	Curved	



Protected by US patent 9,335,455 and other patents and patent applications worldwide

19 GAUGE ILLUMINATION CONNECTOR

A527100	Straight	
A547100	Angled	
A567100	Curved	20g
A587100	Straight, Extended	
A597100	Angled, Extended	
A5A7100	Curved, Extended	
B527103	Straight	
B577103	Angled	23g
B587103	Straight, Extended	
B5B7103	Angled, Extended	
B527104	Straight	
B567104	Curved	25g
B587104	Straight, Extended	
B5A7104	Curved, Extended	
B527105	Straight	27g
B567105	Curved	



B&L STELLARIS PC ILLUMINATION CONNECTOR

A527500	Straight	
A547500	Angled	
A567500	Curved	20g
A587500	Straight, Extended	
A597500	Angled, Extended	
A5A7500	Curved, Extended	
B527503	Straight	
B577503	Angled	23g
B587503	Straight, Extended	
B5B7503	Angled, Extended	
B527504	Straight	
B567504	Curved	25g
B587504	Straight, Extended	
B5A7504	Curved, Extended	
B527505	Straight	27g
B567505	Curved	



ALCON OCUTOME ILLUMINATION CONNECTOR

A527300	Straight	
A547300	Angled	
A567300	Curved	20g
A587300	Straight, Extended	
A597300	Angled, Extended	
A5A7300	Curved, Extended	
B527303	Straight	
B577303	Angled	23g
B587303	Straight, Extended	
B5B7303	Angled, Extended	
B527304	Straight	
B567304	Curved	25g
B587304	Straight, Extended	
B5A7304	Curved, Extended	
B527305	Straight	27g
B567305	Curved	



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

A527400	Straight	
A547400	Angled	
A567400	Curved	20g
A587400	Straight, Extended	
A597400	Angled, Extended	
A5A7400	Curved, Extended	
B527403	Straight	
B577403	Angled	23g
B587403	Straight, Extended	
B5B7403	Angled, Extended	
B527404	Straight	
B567404	Curved	25g
B587404	Straight, Extended	
B5A7404	Curved, Extended	
B527405	Straight	27g
B567405	Curved	



IRIDEX LASER CONNECTOR

ALCON ACCURUS AND CONSTELLATION ILLUMINATION CONNECTOR

A528200	Straight	
A548200	Angled	
A568200	Curved	20g
A588200	Straight, Extended	
A598200	Angled, Extended	
A5A8200	Curved, Extended	
B528203	Straight	
B578203	Angled	23g
B588203	Straight, Extended	
B5B8203	Angled, Extended	
B528204	Straight	
B568204	Curved	25g
B588204	Straight, Extended	
B5A8204	Curved, Extended	
B528205	Straight	27g
B568205	Curved	



Protected by US patent 9,335,455 and other patents and patent applications worldwide

19 GAUGE ILLUMINATION CONNECTOR

A528100	Straight	
A548100	Angled	
A568100	Curved	20g
A588100	Straight, Extended	
A598100	Angled, Extended	
A5A8100	Curved, Extended	
B528103	Straight	
B578103	Angled	23g
B588103	Straight, Extended	
B5B8103	Angled, Extended	
B528104	Straight	
B568104	Curved	25g
B588104	Straight, Extended	
B5A8104	Curved, Extended	
B528105	Straight	27g
B568105	Curved	



B&L STELLARIS PC ILLUMINATION CONNECTOR

A528500	Straight	
A548500	Angled	
A568500	Curved	20g
A588500	Straight, Extended	
A598500	Angled, Extended	
A5A8500	Curved, Extended	
B528503	Straight	
B578503	Angled	23g
B588503	Straight, Extended	
B5B8503	Angled, Extended	
B528504	Straight	
B568504	Curved	25g
B588504	Straight, Extended	
B5A8504	Curved, Extended	
B528505	Straight	27g
B568505	Curved	



ALCON OCUTOME ILLUMINATION CONNECTOR

A528300	Straight	
A548300	Angled	
A568300	Curved	20g
A588300	Straight, Extended	
A598300	Angled, Extended	
A5A8300	Curved, Extended	
B528303	Straight	
B578303	Angled	23g
B588303	Straight, Extended	
B5B8303	Angled, Extended	
B528304	Straight	
B568304	Curved	25g
B588304	Straight, Extended	
B5A8304	Curved, Extended	
B528305	Straight	27g
B568305	Curved	



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

A528400	Straight	
A548400	Angled	
A568400	Curved	20g
A588400	Straight, Extended	
A598400	Angled, Extended	
A5A8400	Curved, Extended	
B528403	Straight	
B578403	Angled	23g
B588403	Straight, Extended	
B5B8403	Angled, Extended	
B528404	Straight	
B568404	Curved	25g
B588404	Straight, Extended	
B5A8404	Curved, Extended	
B528405	Straight	27g
B568405	Curved	



DORC (NON-USA) LASER CONNECTOR

ALCON ACCURUS AND CONSTELLATION ILLUMINATION CONNECTOR

A523200	Straight	
A543200	Angled	
A563200	Curved	20g
A583200	Straight, Extended	
A593200	Angled, Extended	
A5A3200	Curved, Extended	
B523203	Straight	
B573203	Angled	23g
B583203	Straight, Extended	
B5B3203	Angled, Extended	
B523204	Straight	
B563204	Curved	25g
B583204	Straight, Extended	
B5A3204	Curved, Extended	
B523205	Straight	27g
B563205	Curved	



Protected by US patent 9,335,455 and other patents and patent applications worldwide

19 GAUGE ILLUMINATION CONNECTOR

A523100	Straight	
A543100	Angled	
A563100	Curved	20g
A583100	Straight, Extended	
A593100	Angled, Extended	
A5A3100	Curved, Extended	
B523103	Straight	
B573103	Angled	23g
B583103	Straight, Extended	
B5B3103	Angled, Extended	
B523104	Straight	
B563104	Curved	25g
B583104	Straight, Extended	
B5A3104	Curved, Extended	
B523105	Straight	27g
B563105	Curved	



B&L STELLARIS PC ILLUMINATION CONNECTOR

A523500	Straight	
A543500	Angled	
A563500	Curved	20g
A583500	Straight, Extended	
A593500	Angled, Extended	
A5A3500	Curved, Extended	
B523503	Straight	
B573503	Angled	23g
B583503	Straight, Extended	
B5B3503	Angled, Extended	
B523504	Straight	
B563504	Curved	25g
B583504	Straight, Extended	
B5A3504	Curved, Extended	
B523505	Straight	27g
B563505	Curved	



ALCON OCUTOME ILLUMINATION CONNECTOR

A523300	Straight	
A543300	Angled	
A563300	Curved	20g
A583300	Straight, Extended	
A593300	Angled, Extended	
A5A3300	Curved, Extended	
B523303	Straight	
B573303	Angled	23g
B583303	Straight, Extended	
B5B3303	Angled, Extended	
B523304	Straight	
B563304	Curved	25g
B583304	Straight, Extended	
B5A3304	Curved, Extended	
B523305	Straight	27g
B563305	Curved	



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

A523400	Straight	
A543400	Angled	
A563400	Curved	20g
A583400	Straight, Extended	
A593400	Angled, Extended	
A5A3400	Curved, Extended	
B523403	Straight	
B573403	Angled	23g
B583403	Straight, Extended	
B5B3403	Angled, Extended	
B523404	Straight	
B563404	Curved	25g
B583404	Straight, Extended	
B5A3404	Curved, Extended	
B523405	Straight	27g
B563405	Curved	



TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS



QUANTEL VITRA / SUPRA LASER CONNECTOR (DIRECT)

ALCON ACCURUS AND CONSTELLATION ILLUMINATION CONNECTOR

A529200	Straight	
A549200	Angled	
A569200	Curved	20g
A589200	Straight, Extended	
A599200	Angled, Extended	
A5A9200	Curved, Extended	
<hr/>		
B529203	Straight	
B579203	Angled	23g
B589203	Straight, Extended	
B5B9203	Angled, Extended	
<hr/>		
B529204	Straight	
B569204	Curved	25g
B589204	Straight, Extended	
B5A9204	Curved, Extended	
<hr/>		
B529205	Straight	27g
B569205	Curved	



Protected by US patent 9,335,455 and other patents and patent applications worldwide

19 GAUGE ILLUMINATION CONNECTOR

A529100	Straight	
A549100	Angled	
A569100	Curved	20g
A589100	Straight, Extended	
A599100	Angled, Extended	
A5A9100	Curved, Extended	
<hr/>		
B529103	Straight	
B579103	Angled	23g
B589103	Straight, Extended	
B5B9103	Angled, Extended	
<hr/>		
B529104	Straight	
B569104	Curved	25g
B589104	Straight, Extended	
B5A9104	Curved, Extended	
<hr/>		
B529105	Straight	27g
B569105	Curved	



B&L STELLARIS PC ILLUMINATION CONNECTOR

A529500	Straight	
A549500	Angled	
A569500	Curved	20g
A589500	Straight, Extended	
A599500	Angled, Extended	
A5A9500	Curved, Extended	
<hr/>		
B529503	Straight	
B579503	Angled	23g
B589503	Straight, Extended	
B5B9503	Angled, Extended	
<hr/>		
B529504	Straight	
B569504	Curved	25g
B589504	Straight, Extended	
B5A9504	Curved, Extended	
<hr/>		
B529505	Straight	27g
B569505	Curved	



ALCON OCUTOME ILLUMINATION CONNECTOR

A529300	Straight	
A549300	Angled	
A569300	Curved	20g
A589300	Straight, Extended	
A599300	Angled, Extended	
A5A9300	Curved, Extended	
<hr/>		
B529303	Straight	
B579303	Angled	23g
B589303	Straight, Extended	
B5B9303	Angled, Extended	
<hr/>		
B529304	Straight	
B569304	Curved	25g
B589304	Straight, Extended	
B5A9304	Curved, Extended	
<hr/>		
B529305	Straight	27g
B569305	Curved	



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

A529400	Straight	
A549400	Angled	
A569400	Curved	20g
A589400	Straight, Extended	
A599400	Angled, Extended	
A5A9400	Curved, Extended	
<hr/>		
B529403	Straight	
B579403	Angled	23g
B589403	Straight, Extended	
B5B9403	Angled, Extended	
<hr/>		
B529404	Straight	
B569404	Curved	25g
B589404	Straight, Extended	
B5A9404	Curved, Extended	
<hr/>		
B529405	Straight	27g
B569405	Curved	



CE
0413

LUMENIS (NON-HGM) LASER CONNECTOR

ALCON ACCURUS AND CONSTELLATION ILLUMINATION CONNECTOR

A52A200	Straight	
A54A200	Angled	
A56A200	Curved	20g
A58A200	Straight, Extended	
A59A200	Angled, Extended	
A5AA200	Curved, Extended	

B52A203	Straight	
B57A203	Angled	23g
B58A203	Straight, Extended	
B5BA203	Angled, Extended	

B52A204	Straight	
B56A204	Curved	25g
B58A204	Straight, Extended	
B5AA204	Curved, Extended	

B52A205	Straight	27g
B56A205	Curved	



Protected by US patent 9,335,455 and other patents and patent applications worldwide

19 GAUGE ILLUMINATION CONNECTOR

A52A100	Straight	
A54A100	Angled	
A56A100	Curved	20g
A58A100	Straight, Extended	
A59A100	Angled, Extended	
A5AA100	Curved, Extended	

B52A103	Straight	
B57A103	Angled	23g
B58A103	Straight, Extended	
B5BA103	Angled, Extended	

B52A104	Straight	
B56A104	Curved	25g
B58A104	Straight, Extended	
B5AA104	Curved, Extended	

B52A105	Straight	27g
B56A105	Curved	



B&L STELLARIS PC ILLUMINATION CONNECTOR

A52A500	Straight	
A54A500	Angled	
A56A500	Curved	20g
A58A500	Straight, Extended	
A59A500	Angled, Extended	
A5AA500	Curved, Extended	

B52A503	Straight	
B57A503	Angled	23g
B58A503	Straight, Extended	
B5BA503	Angled, Extended	

B52A504	Straight	
B56A504	Curved	25g
B58A504	Straight, Extended	
B5AA504	Curved, Extended	

B52A505	Straight	27g
B56A505	Curved	



ALCON OCUTOME ILLUMINATION CONNECTOR

A52A300	Straight	
A54A300	Angled	
A56A300	Curved	20g
A58A300	Straight, Extended	
A59A300	Angled, Extended	
A5AA300	Curved, Extended	

B52A303	Straight	
B57A303	Angled	23g
B58A303	Straight, Extended	
B5BA303	Angled, Extended	

B52A304	Straight	
B56A304	Curved	25g
B58A304	Straight, Extended	
B5AA304	Curved, Extended	

B52A305	Straight	27g
B56A305	Curved	



A52A400	Straight	
A54A400	Angled	
A56A400	Curved	20g
A58A400	Straight, Extended	
A59A400	Angled, Extended	
A5AA400	Curved, Extended	

B52A403	Straight	
B57A403	Angled	23g
B58A403	Straight, Extended	
B5BA403	Angled, Extended	

B52A404	Straight	
B56A404	Curved	25g
B58A404	Straight, Extended	
B5AA404	Curved, Extended	

B52A405	Straight	27g
B56A405	Curved	



NIDEK LASER CONNECTOR - THREADED

ALCON ACCURUS AND CONSTELLATION ILLUMINATION CONNECTOR

A522200	Straight	
A542200	Angled	
A562200	Curved	20g
A582200	Straight, Extended	
A592200	Angled, Extended	
A5A2200	Curved, Extended	

B522203	Straight	
B572203	Angled	23g
B582203	Straight, Extended	
B5B2203	Angled, Extended	

B522204	Straight	
B562204	Curved	25g
B582204	Straight, Extended	
B5A2204	Curved, Extended	

B522205	Straight	27g
B562205	Curved	



Protected by US patent 9,335,455 and other patents and patent applications worldwide

19 GAUGE ILLUMINATION CONNECTOR

A522100	Straight	
A542100	Angled	
A562100	Curved	20g
A582100	Straight, Extended	
A592100	Angled, Extended	
A5A2100	Curved, Extended	

B522103	Straight	
B572103	Angled	23g
B582103	Straight, Extended	
B5B2103	Angled, Extended	

B522104	Straight	
B562104	Curved	25g
B582104	Straight, Extended	
B5A2104	Curved, Extended	

B522105	Straight	27g
B562105	Curved	



B&L STELLARIS PC ILLUMINATION CONNECTOR

A522500	Straight	
A542500	Angled	
A562500	Curved	20g
A582500	Straight, Extended	
A592500	Angled, Extended	
A5A2500	Curved, Extended	

B522503	Straight	
B572503	Angled	23g
B582503	Straight, Extended	
B5B2503	Angled, Extended	

B522504	Straight	
B562504	Curved	25g
B582504	Straight, Extended	
B5A2504	Curved, Extended	

B522505	Straight	27g
B562505	Curved	



ALCON OCUTOME ILLUMINATION CONNECTOR

A522300	Straight	
A542300	Angled	
A562300	Curved	20g
A582300	Straight, Extended	
A592300	Angled, Extended	
A5A2300	Curved, Extended	

B522303	Straight	
B572303	Angled	23g
B582303	Straight, Extended	
B5B2303	Angled, Extended	

B522304	Straight	
B562304	Curved	25g
B582304	Straight, Extended	
B5A2304	Curved, Extended	

B522305	Straight	27g
B562305	Curved	



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

A522400	Straight	
A542400	Angled	
A562400	Curved	20g
A582400	Straight, Extended	
A592400	Angled, Extended	
A5A2400	Curved, Extended	

B522403	Straight	
B572403	Angled	23g
B582403	Straight, Extended	
B5B2403	Angled, Extended	

B522404	Straight	
B562404	Curved	25g
B582404	Straight, Extended	
B5A2404	Curved, Extended	

B522405	Straight	27g
B562405	Curved	



NIDEK LASER CONNECTOR - PUSH-IN

ALCON ACCURUS AND CONSTELLATION ILLUMINATION CONNECTOR

A525200	Straight	
A545200	Angled	
A565200	Curved	20g
A585200	Straight, Extended	
A595200	Angled, Extended	
A5A5200	Curved, Extended	
<hr/>		
B525203	Straight	
B575203	Angled	23g
B585203	Straight, Extended	
B5B5203	Angled, Extended	
<hr/>		
B525204	Straight	
B565204	Curved	25g
B585204	Straight, Extended	
B5A5204	Curved, Extended	



Protected by US patent 9,335,455 and other patents and patent applications worldwide

19 GAUGE ILLUMINATION CONNECTOR

A525100	Straight	
A545100	Angled	
A565100	Curved	20g
A585100	Straight, Extended	
A595100	Angled, Extended	
A5A5100	Curved, Extended	
<hr/>		
B525103	Straight	
B575103	Angled	23g
B585103	Straight, Extended	
B5B5103	Angled, Extended	
<hr/>		
B525104	Straight	
B565104	Curved	25g
B585104	Straight, Extended	
B5A5104	Curved, Extended	



B&L STELLARIS PC ILLUMINATION CONNECTOR

A525500	Straight	
A545500	Angled	
A565500	Curved	20g
A585500	Straight, Extended	
A595500	Angled, Extended	
A5A5500	Curved, Extended	
<hr/>		
B525503	Straight	
B575503	Angled	23g
B585503	Straight, Extended	
B5B5503	Angled, Extended	
<hr/>		
B525504	Straight	
B565504	Curved	25g
B585504	Straight, Extended	
B5A5504	Curved, Extended	



ALCON OCUTOME ILLUMINATION CONNECTOR

A525300	Straight	
A545300	Angled	
A565300	Curved	20g
A585300	Straight, Extended	
A595300	Angled, Extended	
A5A5300	Curved, Extended	
<hr/>		
B525303	Straight	
B575303	Angled	23g
B585303	Straight, Extended	
B5B5303	Angled, Extended	
<hr/>		
B525304	Straight	
B565304	Curved	25g
B585304	Straight, Extended	
B5A5304	Curved, Extended	



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

A525400	Straight	
A545400	Angled	
A565400	Curved	20g
A585400	Straight, Extended	
A595400	Angled, Extended	
A5A5400	Curved, Extended	
<hr/>		
B525403	Straight	
B575403	Angled	23g
B585403	Straight, Extended	
B5B5403	Angled, Extended	
<hr/>		
B525404	Straight	
B565404	Curved	25g
B585404	Straight, Extended	
B5A5404	Curved, Extended	



LIGHTMED LASER CONNECTOR

ALCON ACCURUS AND CONSTELLATION ILLUMINATION CONNECTOR

A52B200	Straight	
A54B200	Angled	
A56B200	Curved	20g
A58B200	Straight, Extended	
A59B200	Angled, Extended	
A5AB200	Curved, Extended	

B52B203	Straight	
B57B203	Angled	23g
B58B203	Straight, Extended	
B5BB203	Angled, Extended	

B52B204	Straight	
B56B204	Curved	25g
B58B204	Straight, Extended	
B5AB204	Curved, Extended	

B52B205	Straight	27g
B56B205	Curved	



Protected by US patent 9,335,455 and other patents and patent applications worldwide

19 GAUGE ILLUMINATION CONNECTOR

A52B100	Straight	
A54B100	Angled	
A56B100	Curved	20g
A58B100	Straight, Extended	
A59B100	Angled, Extended	
A5AB100	Curved, Extended	

B52B103	Straight	
B57B103	Angled	23g
B58B103	Straight, Extended	
B5BB103	Angled, Extended	

B52B104	Straight	
B56B104	Curved	25g
B58B104	Straight, Extended	
B5AB104	Curved, Extended	

B52B105	Straight	27g
B56B105	Curved	



B&L STELLARIS PC ILLUMINATION CONNECTOR

A52B500	Straight	
A54B500	Angled	
A56B500	Curved	20g
A58B500	Straight, Extended	
A59B500	Angled, Extended	
A5AB500	Curved, Extended	

B52B503	Straight	
B57B503	Angled	23g
B58B503	Straight, Extended	
B5BB503	Angled, Extended	

B52B504	Straight	
B56B504	Curved	25g
B58B504	Straight, Extended	
B5AB504	Curved, Extended	

B52B505	Straight	27g
B56B505	Curved	



ALCON OCUTOME ILLUMINATION CONNECTOR

A52B300	Straight	
A54B300	Angled	
A56B300	Curved	20g
A58B300	Straight, Extended	
A59B300	Angled, Extended	
A5AB300	Curved, Extended	

B52B303	Straight	
B57B303	Angled	23g
B58B303	Straight, Extended	
B5BB303	Angled, Extended	

B52B304	Straight	
B56B304	Curved	25g
B58B304	Straight, Extended	
B5AB304	Curved, Extended	

B52B305	Straight	27g
B56B305	Curved	



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

A52B400	Straight	
A54B400	Angled	
A56B400	Curved	20g
A58B400	Straight, Extended	
A59B400	Angled, Extended	
A5AB400	Curved, Extended	

B52B403	Straight	
B57B403	Angled	23g
B58B403	Straight, Extended	
B5BB403	Angled, Extended	

B52B404	Straight	
B56B404	Curved	25g
B58B404	Straight, Extended	
B5AB404	Curved, Extended	

B52B405	Straight	27g
B56B405	Curved	



HGM LASER CONNECTOR - NO BAR CODE

ALCON ACCURUS AND CONSTELLATION ILLUMINATION CONNECTOR

A524200	Straight	
A544200	Angled	
A564200	Curved	20g
A584200	Straight, Extended	
A594200	Angled, Extended	
A5A4200	Curved, Extended	
B524203	Straight	
B574203	Angled	23g
B584203	Straight, Extended	
B5B4203	Angled, Extended	
B524204	Straight	
B564204	Curved	25g
B584204	Straight, Extended	
B5A4204	Curved, Extended	



Protected by US patent 9,335,455 and other patents and patent applications worldwide

19 GAUGE ILLUMINATION CONNECTOR

A524100	Straight	
A544100	Angled	
A564100	Curved	20g
A584100	Straight, Extended	
A594100	Angled, Extended	
A5A4100	Curved, Extended	
B524103	Straight	
B574103	Angled	23g
B584103	Straight, Extended	
B5B4103	Angled, Extended	
B524104	Straight	
B564104	Curved	25g
B584104	Straight, Extended	
B5A4104	Curved, Extended	



B&L STELLARIS PC ILLUMINATION CONNECTOR

A524500	Straight	
A544500	Angled	
A564500	Curved	20g
A584500	Straight, Extended	
A594500	Angled, Extended	
A5A4500	Curved, Extended	
B524503	Straight	
B574503	Angled	23g
B584503	Straight, Extended	
B5B4503	Angled, Extended	
B524504	Straight	
B564504	Curved	25g
B584504	Straight, Extended	
B5A4504	Curved, Extended	



ALCON OCUTOME ILLUMINATION CONNECTOR

A524300	Straight	
A544300	Angled	
A564300	Curved	20g
A584300	Straight, Extended	
A594300	Angled, Extended	
A5A4300	Curved, Extended	
B524303	Straight	
B574303	Angled	23g
B584303	Straight, Extended	
B5B4303	Angled, Extended	
B524304	Straight	
B564304	Curved	25g
B584304	Straight, Extended	
B5A4304	Curved, Extended	



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

A524400	Straight	
A544400	Angled	
A564400	Curved	20g
A584400	Straight, Extended	
A594400	Angled, Extended	
A5A4400	Curved, Extended	
B524403	Straight	
B574403	Angled	23g
B584403	Straight, Extended	
B5B4403	Angled, Extended	
B524404	Straight	
B564404	Curved	25g
B584404	Straight, Extended	
B5A4404	Curved, Extended	





Aspiration add-on functionality is available upon request.

A flexible plastic jacket covers and protects the 8ft long fiber.

Homogeneous energy distribution as well as round and consistently small laser spot allow the tip to stay few millimeters away from the retina.

Ergonomic and non-slip handpiece design.



Different connectors that suit most laser and illumination machines.

Sterilized by Ethylene Oxide, double pouched, single-use, latex-free.

Available in Dia 0.41 mm (27ga), Dia 0.51 mm (25ga), Dia 0.64 mm (23ga), and Dia 0.91 mm (20g). Length 2.4 m (8 ft). Box of 10.

Do not re-sterilize. Do not reuse.

Material fatigue decreases performance.

WARNING

Follow the recommendations of the laser and illumination source manufacturers (e.g.: calibration of the laser source as well as the need for protective eye wear for the corresponding wavelength of the laser model used).

WARNING IN CASE OF ILLUMINATION

Due to the potentially hazardous effects of phototoxicity to the retina from prolonged endoillumination, this device should be used with an illumination source with filters that eliminate UV radiation (<400nm) and whenever possible, with filters that eliminate short wavelength blue light (<420nm). Care should be taken not to focus this device on a single point of the retina for prolonged, unnecessary periods of time.

CAUTION: FEDERAL (U.S.A) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN

CONTACT INFORMATION:

+1 770 777 6613

+1 877 4 25 GAUGE (4284)

+1 678 623 3765 (FAX)

info@ophthalmed.com

1050 Northfield Court, Suite 280
Roswell, GA 30076 USA

DISTRIBUTED BY:



OPHTHALMEDTM

MADE IN THE USA