




OPHTHALMEDTM

LASER PROBES
FOR ENDO-OCULAR
PHOTOCOAGULATION
TREATMENTS

WWW.OPHTHALMED.COM

+1 877 4 25 GAUGE (4284)

In business since 2001, Ophthalmed LLC specializes in designing and manufacturing higher quality endo photocoagulation laser probes.

The range of our products is one of the most extensive on the market, backed up by US and European patents, and private labeled to major companies in US and worldwide, including well known laser machine manufacturers and other high profile distributors.

We are committed to serving the needs of the ophthalmic industry, and to providing the latest technologies to the surgeons, including custom designs.

Our current product range includes laser and illuminated laser probes in various tip sizes and configurations, including bending (patented). We connect to most laser and illumination machines.

Our probes are manufactured in accordance with a quality system certified to ISO 13485 and cleared/approved by various regulatory agencies, such as US FDA and the European Community.

PRODUCT GUIDE

PAGE

Laser probes

1-2

Straight, angled and curved tips, 20g to 27g, 9 laser connectors

Bending laser probes

3

90° bending tip, 20g to 27g, 7 laser connectors

Illuminated laser probes

4-12

Straight, angled and curved tips, 20g to 27g, 9 laser connectors, 6 illumination connectors

Illuminated bending laser probes

13-16

90° bending tip, 20g to 27g, 7 laser connectors, 6 illumination connectors

Product index

17-20

All 714 codes in code order, with their page numbers

	STRAIGHT	ANGLED	CURVED
20 g			
23 g			
25g			
27g			

ALCON / ZEISS CONNECTOR

A217000	Straight	
A237000	Angled	20g
A257000	Curved	
A827000	Straight	23g
A877000	Angled	
A327000	Straight	
A377000	Angled	25g
A367000	Curved	
AA27000	Straight	27g
AA67000	Curved	

**IRIDEX CONNECTOR**

A218000	Straight	
A238000	Angled	20g
A258000	Curved	
A828000	Straight	23g
A878000	Angled	
A328000	Straight	
A378000	Angled	25g
A368000	Curved	
AA28000	Straight	27g
AA68000	Curved	

**DORC (NON-USA) CONNECTOR**

A213000	Straight	
A233000	Angled	20g
A253000	Curved	
A823000	Straight	23g
A873000	Angled	
A323000	Straight	
A373000*	Angled	25g
A363000	Curved	
AA23000	Straight	27g
AA63000	Curved	

**QUANTEL SUPRA/VITRA CONNECTOR (DIRECT)**

A219000	Straight	
A239000	Angled	20g
A259000	Curved	
A829000	Straight	23g
A879000	Angled	
A329000	Straight	
A379000	Angled	25g
A369000	Curved	
AA29000	Straight	27g
AA69000	Curved	

**LUMENIS (NON-HGM) CONNECTOR**

A21A000	Straight	
A23A000	Angled	20g
A25A000	Curved	
A82A000	Straight	23g
A87A000	Angled	
A32A000	Straight	
A37A000	Angled	25g
A36A000	Curved	
AA2A000	Straight	27g
AA6A000	Curved	

**NIDEK CONNECTOR - THREADED**

A212000	Straight	
A232000	Angled	20g
A252000	Curved	
A822000	Straight	23g
A872000	Angled	
A322000	Straight	
A372000	Angled	25g
A362000	Curved	
AA22000	Straight	27g
AA62000	Curved	



ADDITIONAL VARIATIONS ARE AVAILABLE.

* Not CE marked yet

CE
0413

TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS

	STRAIGHT	ANGLED	CURVED
20 g			
23 g			
25g			
27g			

NIDEK CONNECTOR - PUSH-IN

A215000	Straight	
A235000	Angled	20g
A255000	Curved	
A825000	Straight	23g
A875000	Angled	
A325000	Straight	
A375000	Angled	25g
A365000	Curved	
AA25000	Straight	27g
AA65000	Curved	



LIGHTMED CONNECTOR

A21B000	Straight	
A23B000	Angled	20g
A25B000	Curved	
A82B000	Straight	23g
A87B000	Angled	
A32B000	Straight	
A37B000	Angled	25g
A36B000	Curved	
AA2B000	Straight	27g
AA6B000	Curved	



HGM CONNECTOR - NO BAR CODE


A214000	Straight	
A234000	Angled	20g
A254000	Curved	
A824000	Straight	23g
A874000	Angled	
A324000	Straight	
A374000	Angled	25g
A364000	Curved	
AA24000	Straight	27g
AA64000	Curved	



BENDING LASER PROBES

Protected by US Pat. 9,370,447, EP 2765943
and corresponding national patents.


Bending probe features: see page 13.



Fully slide button forward to straighten the tip prior to penetration or withdrawal; sliding the button all the way forward straightens the tip while sliding the button backward incrementally curves the tip.

REF	TIP	LASER CONNECTOR	
A207000	BENDING, 20G	ALCON / ZEISS	
A807000	BENDING, 23G		
A307000	BENDING, 25G		
AA07000	BENDING, 27G		
A208000	BENDING, 20G	IRIDEX	
A808000	BENDING, 23G		
A308000	BENDING, 25G		
AA08000	BENDING, 27G		
A203000	BENDING, 20G	DORC (non-USA)	
A803000	BENDING, 23G		
A303000	BENDING, 25G		
AA03000	BENDING, 27G		
A209000	BENDING, 20G	QUANTEL SUPRA/VITRA (DIRECT CONNECTION)	
A809000	BENDING, 23G		
A309000	BENDING, 25G		
AA09000	BENDING, 27G		
A20A000	BENDING, 20G	LUMENIS (non-HGM)	
A80A000	BENDING, 23G		
A30A000	BENDING, 25G		
AA0A000	BENDING, 27G		
A202000	BENDING, 20G	NIDEK (THREADED)	
A802000	BENDING, 23G		
A302000	BENDING, 25G		
AA02000	BENDING, 27G		
A20B000	BENDING, 20G	LIGHTMED	
A80B000	BENDING, 23G		
A30B000	BENDING, 25G		
AA0B000	BENDING, 27G		

STANDARD TIPS



ALCON / ZEISS LASER CONNECTOR

ALCON STANDARD ILLUMINATION CONNECTOR

A527200	Straight	20g
A547200	Angled	
A567200	Curved	
B527203	Straight	23g
B577203	Angled	
B527204	Straight	25g
B577204	Angled	
B567204	Curved	
B527205	Straight	27g
B567205	Curved	



B&L STELLARIS PC ILLUMINATION CONNECTOR

A527500	Straight	20g
A547500	Angled	
A567500	Curved	
B527503	Straight	23g
B577503	Angled	
B527504	Straight	25g
B577504	Angled	
B567504	Curved	
B527505	Straight	27g
B567505	Curved	



19 GAUGE ILLUMINATION CONNECTOR

A527100	Straight	20g
A547100	Angled	
A567100	Curved	
B527103	Straight	23g
B577103	Angled	
B527104	Straight	25g
B577104	Angled	
B567104	Curved	
B527105	Straight	27g
B567105	Curved	



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

A527400	Straight	20g
A547400	Angled	
A567400	Curved	
B527403	Straight	23g
B577403	Angled	
B527404	Straight	25g
B577404	Angled	
B567404	Curved	
B527405	Straight	27g
B567405	Curved	



ALCON OCUTOME ILLUMINATION CONNECTOR

A527300	Straight	20g
A547300	Angled	
A567300	Curved	
B527303	Straight	23g
B577303	Angled	
B527304	Straight	25g
B577304	Angled	
B567304	Curved	
B527305	Straight	27g
B567305	Curved	



DORC EVA ILLUMINATION CONNECTOR

A527600	Straight	20g
A547600	Angled	
A567600	Curved	
B527603	Straight	23g
B577603	Angled	
B527604	Straight	25g
B577604	Angled	
B567604	Curved	
B527605	Straight	27g
B567605	Curved	



CE
0413

TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS

IRIDEX LASER CONNECTOR

**ALCON STANDARD
ILLUMINATION CONNECTOR**

A528200	Straight	
A548200	Angled	20g
A568200	Curved	
B528203	Straight	
B578203	Angled	23g
B528204	Straight	
B578204	Angled	25g
B568204	Curved	
B528205	Straight	
B568205	Curved	27g

**19 GAUGE
ILLUMINATION CONNECTOR**

A528100	Straight	
A548100	Angled	20g
A568100	Curved	
B528103	Straight	
B578103	Angled	23g
B528104	Straight	
B578104	Angled	25g
B568104	Curved	
B528105	Straight	
B568105	Curved	27g

**ALCON OCUTOME
ILLUMINATION CONNECTOR**

A528300	Straight	
A548300	Angled	20g
A568300	Curved	
B528303	Straight	
B578303	Angled	23g
B528304	Straight	
B578304	Angled	25g
B568304	Curved	
B528305	Straight	
B568305	Curved	27g

**B&L STELLARIS PC
ILLUMINATION CONNECTOR**

A528500	Straight	
A548500	Angled	20g
A568500	Curved	
B528503	Straight	
B578503	Angled	23g
B528504	Straight	
B578504	Angled	25g
B568504	Curved	
B528505	Straight	
B568505	Curved	27g

**B&L / STORZ MILLENIUM
ILLUMINATION CONNECTOR**

A528400	Straight	
A548400	Angled	20g
A568400	Curved	
B528403	Straight	
B578403	Angled	23g
B528404	Straight	
B578404	Angled	25g
B568404	Curved	
B528405	Straight	
B568405	Curved	27g

**DORC EVA
ILLUMINATION CONNECTOR**

A528600	Straight	
A548600	Angled	20g
A568600	Curved	
B528603	Straight	
B578603	Angled	23g
B528604	Straight	
B578604	Angled	25g
B568604	Curved	
B528605	Straight	
B568605	Curved	27g



0413

DORC (NON-USA) LASER CONNECTOR

ALCON STANDARD ILLUMINATION CONNECTOR

A523200	Straight	20g
A543200	Angled	
A563200	Curved	
B523203	Straight	23g
B573203	Angled	
B523204	Straight	25g
B573204	Angled	
B563204	Curved	
B523205	Straight	27g
B563205	Curved	



B&L STELLARIS PC ILLUMINATION CONNECTOR

A523500	Straight	20g
A543500	Angled	
A563500	Curved	
B523503	Straight	23g
B573503	Angled	
B523504	Straight	25g
B573504	Angled	
B563504	Curved	
B523505	Straight	27g
B563505	Curved	



19 GAUGE ILLUMINATION CONNECTOR

A523100	Straight	20g
A543100	Angled	
A563100	Curved	
B523103	Straight	23g
B573103	Angled	
B523104	Straight	25g
B573104	Angled	
B563104	Curved	
B523105	Straight	27g
B563105	Curved	



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

A523400	Straight	20g
A543400	Angled	
A563400	Curved	
B523403	Straight	23g
B573403	Angled	
B523404	Straight	25g
B573404	Angled	
B563404	Curved	
B523405	Straight	27g
B563405	Curved	



ALCON OCUTOME ILLUMINATION CONNECTOR

A523300	Straight	20g
A543300	Angled	
A563300	Curved	
B523303	Straight	23g
B573303	Angled	
B523304	Straight	25g
B573304	Angled	
B563304	Curved	
B523305	Straight	27g
B563305	Curved	



DORC EVA ILLUMINATION CONNECTOR

A523600	Straight	20g
A543600	Angled	
A563600	Curved	
B523603	Straight	23g
B573603	Angled	
B523604	Straight	25g
B573604	Angled	
B563604	Curved	
B523605	Straight	27g
B563605	Curved	



CE
0413

TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS

QUANTEL SUPRA/VITRA LASER CONNECTOR (DIRECT)

ALCON STANDARD
ILLUMINATION CONNECTOR

A529200	Straight	
A549200	Angled	20g
A569200	Curved	
B529203	Straight	
B579203	Angled	23g
B529204	Straight	
B579204	Angled	25g
B569204	Curved	
B529205	Straight	
B569205	Curved	27g

19 GAUGE
ILLUMINATION CONNECTOR

A529100	Straight	
A549100	Angled	20g
A569100	Curved	
B529103	Straight	
B579103	Angled	23g
B529104	Straight	
B579104	Angled	25g
B569104	Curved	
B529105	Straight	
B569105	Curved	27g

ALCON OCUTOME
ILLUMINATION CONNECTOR

A529300	Straight	
A549300	Angled	20g
A569300	Curved	
B529303	Straight	
B579303	Angled	23g
B529304	Straight	
B579304	Angled	25g
B569304	Curved	
B529305	Straight	
B569305	Curved	27g

B&L STELLARIS PC
ILLUMINATION CONNECTOR

A529500	Straight	
A549500	Angled	20g
A569500	Curved	
B529503	Straight	
B579503	Angled	23g
B529504	Straight	
B579504	Angled	25g
B569504	Curved	
B529505	Straight	
B569505	Curved	27g

B&L / STORZ MILLENIUM
ILLUMINATION CONNECTOR

A529400	Straight	
A549400	Angled	20g
A569400	Curved	
B529403	Straight	
B579403	Angled	23g
B529404	Straight	
B579404	Angled	25g
B569404	Curved	
B529405	Straight	
B569405	Curved	27g

DORC EVA
ILLUMINATION CONNECTOR

A529600	Straight	
A549600	Angled	20g
A569600	Curved	
B529603	Straight	
B579603	Angled	23g
B529604	Straight	
B579604	Angled	25g
B569604	Curved	
B529605	Straight	
B569605	Curved	27g



LUMENIS (NON-HGM) LASER CONNECTOR

ALCON STANDARD ILLUMINATION CONNECTOR

A52A200	Straight	20g
A54A200	Angled	
A56A200	Curved	
B52A203	Straight	23g
B57A203	Angled	
B52A204	Straight	25g
B57A204	Angled	
B56A204	Curved	
B52A205	Straight	27g
B56A205	Curved	



B&L STELLARIS PC ILLUMINATION CONNECTOR

A52A500	Straight	20g
A54A500	Angled	
A56A500	Curved	
B52A503	Straight	23g
B57A503	Angled	
B52A504	Straight	25g
B57A504	Angled	
B56A504	Curved	
B52A505	Straight	27g
B56A505	Curved	



19 GAUGE ILLUMINATION CONNECTOR

A52A100	Straight	20g
A54A100	Angled	
A56A100	Curved	
B52A103	Straight	23g
B57A103	Angled	
B52A104	Straight	25g
B57A104	Angled	
B56A104	Curved	
B52A105	Straight	27g
B56A105	Curved	



ALCON OCUTOME ILLUMINATION CONNECTOR

A52A300	Straight	20g
A54A300	Angled	
A56A300	Curved	
B52A303	Straight	23g
B57A303*	Angled	
B52A304	Straight	25g
B57A304	Angled	
B56A304	Curved	
B52A305	Straight	27g
B56A305	Curved	



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

A52A400	Straight	20g
A54A400	Angled	
A56A400	Curved	
B52A403	Straight	23g
B57A403	Angled	
B52A404	Straight	25g
B57A404	Angled	
B56A404	Curved	
B52A405	Straight	27g
B56A405	Curved	



DORC EVA ILLUMINATION CONNECTOR

A52A600	Straight	20g
A54A600	Angled	
A56A600	Curved	
B52A603	Straight	23g
B57A603	Angled	
B52A604	Straight	25g
B57A604	Angled	
B56A604	Curved	
B52A605	Straight	27g
B56A605	Curved	



* Not CE marked yet

NIDEK LASER CONNECTOR - THREADED

ALCON STANDARD
ILLUMINATION CONNECTOR

A522200	Straight	
A542200	Angled	20g
A562200	Curved	
B522203	Straight	
B572203	Angled	23g
B522204	Straight	
B572204	Angled	25g
B562204	Curved	
B522205	Straight	
B562205	Curved	27g

19 GAUGE
ILLUMINATION CONNECTOR

A522100	Straight	
A542100	Angled	20g
A562100	Curved	
B522103	Straight	
B572103	Angled	23g
B522104	Straight	
B572104	Angled	25g
B562104	Curved	
B522105	Straight	
B562105	Curved	27g

ALCON OCUTOME
ILLUMINATION CONNECTOR

A522300	Straight	
A542300	Angled	20g
A562300	Curved	
B522303	Straight	
B572303	Angled	23g
B522304	Straight	
B572304	Angled	25g
B562304	Curved	
B522305	Straight	
B562305	Curved	27g

B&L STELLARIS PC
ILLUMINATION CONNECTOR

A522500	Straight	
A542500	Angled	20g
A562500	Curved	
B522503	Straight	
B572503	Angled	23g
B522504	Straight	
B572504	Angled	25g
B562504	Curved	
B522505	Straight	
B562505	Curved	27g

B&L / STORZ MILLENIUM
ILLUMINATION CONNECTOR

A522400	Straight	
A542400	Angled	20g
A562400	Curved	
B522403	Straight	
B572403	Angled	23g
B522404	Straight	
B572404	Angled	25g
B562404	Curved	
B522405	Straight	
B562405	Curved	27g

DORC EVA
ILLUMINATION CONNECTOR

A522600	Straight	
A542600	Angled	20g
A562600	Curved	
B522603	Straight	
B572603	Angled	23g
B522604	Straight	
B572604	Angled	25g
B562604	Curved	
B522605	Straight	
B562605	Curved	27g



LIGHTMED LASER CONNECTOR

ALCON STANDARD ILLUMINATION CONNECTOR

A52B200	Straight	20g
A54B200	Angled	
A56B200	Curved	
B52B203	Straight	23g
B57B203	Angled	
B52B204	Straight	25g
B57B204	Angled	
B56B204	Curved	
B52B205	Straight	27g
B56B205	Curved	



B&L STELLARIS PC ILLUMINATION CONNECTOR

A52B500	Straight	20g
A54B500	Angled	
A56B500	Curved	
B52B503	Straight	23g
B57B503	Angled	
B52B504	Straight	25g
B57B504	Angled	
B56B504	Curved	
B52B505	Straight	27g
B56B505	Curved	



19 GAUGE ILLUMINATION CONNECTOR

A52B100	Straight	20g
A54B100	Angled	
A56B100	Curved	
B52B103	Straight	23g
B57B103	Angled	
B52B104	Straight	25g
B57B104	Angled	
B56B104	Curved	
B52B105	Straight	27g
B56B105	Curved	



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

A52B400	Straight	20g
A54B400	Angled	
A56B400	Curved	
B52B403	Straight	23g
B57B403	Angled	
B52B404	Straight	25g
B57B404*	Angled	
B56B404	Curved	
B52B405	Straight	27g
B56B405	Curved	



ALCON OCUTOME ILLUMINATION CONNECTOR

A52B300	Straight	20g
A54B300	Angled	
A56B300	Curved	
B52B303	Straight	23g
B57B303*	Angled	
B52B304	Straight	25g
B57B304	Angled	
B56B304	Curved	
B52B305	Straight	27g
B56B305	Curved	



DORC EVA ILLUMINATION CONNECTOR

A52B600	Straight	20g
A54B600	Angled	
A56B600	Curved	
B52B603	Straight	23g
B57B603	Angled	
B52B604	Straight	25g
B57B604	Angled	
B56B604	Curved	
B52B605	Straight	27g
B56B605	Curved	



CE
0413

* Not CE marked yet

TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS

HGM LASER CONNECTOR - NO BAR CODE

ALCON STANDARD
ILLUMINATION CONNECTOR

A524200	Straight	
A544200	Angled	20g
A564200	Curved	
B524203	Straight	23g

19 GAUGE
ILLUMINATION CONNECTOR

A524100	Straight	
A544100	Angled	20g
A564100	Curved	
B524103	Straight	23g

ALCON OCUTOME
ILLUMINATION CONNECTOR

A524300	Straight	
A544300	Angled	20g
A564300	Curved	
B524303	Straight	23g

B&L STELLARIS PC
ILLUMINATION CONNECTOR

A524500	Straight	
A544500	Angled	20g
A564500	Curved	
B524503	Straight	23g

B&L / STORZ MILLENIUM
ILLUMINATION CONNECTOR

A524400	Straight	
A544400	Angled	20g
A564400	Curved	
B524403	Straight	23g

DORC EVA
ILLUMINATION CONNECTOR

A524600	Straight	
A544600	Angled	20g
A564600	Curved	
B524603	Straight	23g



NIDEK LASER CONNECTOR - PUSH-IN

ALCON STANDARD
ILLUMINATION CONNECTOR

A525200	Straight	
A545200	Angled	20g
A565200	Curved	
B525203	Straight	23g

19 GAUGE
ILLUMINATION CONNECTOR

A525100	Straight	
A545100	Angled	20g
A565100	Curved	
B525103	Straight	23g

ALCON OCUTOME
ILLUMINATION CONNECTOR

A525300	Straight	
A545300	Angled	20g
A565300	Curved	
B525303	Straight	23g

B&L STELLARIS PC
ILLUMINATION CONNECTOR

A525500	Straight	
A545500	Angled	20g
A565500	Curved	
B525503	Straight	23g

B&L / STORZ MILLENIUM
ILLUMINATION CONNECTOR

A525400	Straight	
A545400	Angled	20g
A565400	Curved	
B525403	Straight	23g

DORC EVA
ILLUMINATION CONNECTOR

A525600	Straight	
A545600	Angled	20g
A565600	Curved	
B525603	Straight	23g




BENDING LASER PROBES



Patented† 90° bending laser probes down to 27g, allowing the surgeon to selectively change tip angle during surgery. Illuminated bending versions with bright widefield illumination are available down to 25g, and in 27g with the Alcon / Zeiss laser connector: B507205 (Alcon Standard illumination) and B507305 (Alcon Ocutome illumination).

FEATURES



- Adjusting angle direction between 0 to 90 degrees
- No externally moving parts. Straightening tube extends from within pre-curved tip
- Stiffer shaft for better maneuverability
- Very easy penetration and withdrawal of flexible tip
- No need to squeeze handle to maintain angle. Tip remains in position without having to hold the sliding button
- No abrupt movement of tip. Angle movement is linear with respect to sliding button
- No risk from tip advancing when changing the angle. Tip remains fixed with respect to surgical site
- **Bright widefield illumination exceeding 80°**
- Extremely uniform energy distribution throughout a small laser spot
- Connect to most laser and illumination machines

† Protected by US Pat. 9,370,447, EP 2765943 and corresponding national patents.

ALCON / ZEISS LASER CONNECTOR



Protected by US Pat. 9,370,447, EP 2765943 and corresponding national patents.

ALCON STANDARD
ILLUMINATION CONNECTOR

A507200	20 g
B507203	23 g
B507204	25 g
B507205	27 g

19 GAUGE
ILLUMINATION CONNECTOR

A507100	20 g
B507103	23 g
B507104	25 g

ALCON OCUTOME
ILLUMINATION CONNECTOR

A507300	20 g
B507303	23 g
B507304	25 g
B507305	27 g

B&L STELLARIS PC
ILLUMINATION CONNECTOR

A507500	20 g
B507503	23 g
B507504	25 g

B&L / STORZ MILLENIUM
ILLUMINATION CONNECTOR

A507400	20 g
B507403	23 g
B507404	25 g

DORC EVA
ILLUMINATION CONNECTOR

A507600	20 g
B507603	23 g
B507604	25 g



IRIDEX LASER CONNECTOR



Protected by US Pat. 9,370,447, EP 2765943 and corresponding national patents.

ALCON STANDARD
ILLUMINATION CONNECTOR

A508200	20 g
B508203	23 g
B508204	25 g

19 GAUGE
ILLUMINATION CONNECTOR

A508100	20 g
B508103	23 g
B508104	25 g

ALCON OCUTOME
ILLUMINATION CONNECTOR

A508300	20 g
B508303	23 g
B508304	25 g

B&L STELLARIS PC
ILLUMINATION CONNECTOR

A508500	20 g
B508503	23 g
B508504	25 g

B&L / STORZ MILLENIUM
ILLUMINATION CONNECTOR

A508400	20 g
B508403	23 g
B508404	25 g

DORC EVA
ILLUMINATION CONNECTOR

A508600	20 g
B508603	23 g
B508604	25 g



DORC (NON-USA) LASER CONNECTOR



Protected by US Pat. 9,370,447, EP 2765943 and corresponding national patents.

ALCON STANDARD ILLUMINATION CONNECTOR

A503200 20 g
B503203 23 g
B503204 25 g



19 GAUGE ILLUMINATION CONNECTOR

A503100 20 g
B503103 23 g
B503104 25 g



ALCON OCUTOME ILLUMINATION CONNECTOR

A503300 20 g
B503303 23 g
B503304 25 g



B&L STELLARIS PC ILLUMINATION CONNECTOR

A503500 20 g
B503503 23 g
B503504 25 g



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

A503400 20 g
B503403 23 g
B503404 25 g



DORC EVA ILLUMINATION CONNECTOR

A503600 20 g
B503603 23 g
B503604 25 g



QUANTEL SUPRA/VITRA LASER CONNECTOR (DIRECT)



Protected by US Pat. 9,370,447, EP 2765943 and corresponding national patents.

ALCON STANDARD ILLUMINATION CONNECTOR

A509200 20 g
B509203 23 g
B509204 25 g



19 GAUGE ILLUMINATION CONNECTOR

A509100 20 g
B509103 23 g
B509104 25 g



ALCON OCUTOME ILLUMINATION CONNECTOR

A509300 20 g
B509303 23 g
B509304 25 g



B&L STELLARIS PC ILLUMINATION CONNECTOR

A509500 20 g
B509503 23 g
B509504 25 g



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

A509400 20 g
B509403 23 g
B509404 25 g



DORC EVA ILLUMINATION CONNECTOR

A509600 20 g
B509603 23 g
B509604 25 g



LUMENIS (NON-HGM) LASER CONNECTOR



Protected by US Pat. 9,370,447, EP 2765943 and corresponding national patents.

ALCON STANDARD ILLUMINATION CONNECTOR

A50A200 20 g
B50A203 23 g
B50A204 25 g



19 GAUGE ILLUMINATION CONNECTOR

A50A100 20 g
B50A103 23 g
B50A104 25 g



ALCON OCUTOME ILLUMINATION CONNECTOR

A50A300 20 g
B50A303 23 g
B50A304 25 g



B&L STELLARIS PC ILLUMINATION CONNECTOR

A50A500 20 g
B50A503 23 g
B50A504 25 g



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

A50A400 20 g
B50A403 23 g
B50A404 25 g



DORC EVA ILLUMINATION CONNECTOR

A50A600 20 g
B50A603 23 g
B50A604 25 g



CE
0413

NIDEK LASER CONNECTOR - THREADED



Protected by US Pat. 9,370,447, EP 2765943 and corresponding national patents.

ALCON STANDARD ILLUMINATION CONNECTOR

A502200	20 g
B502203	23 g
B502204	25 g



19 GAUGE ILLUMINATION CONNECTOR

A502100	20 g
B502103	23 g
B502104	25 g



ALCON OCUTOME ILLUMINATION CONNECTOR

A502300	20 g
B502303	23 g
B502304	25 g



B&L STELLARIS PC ILLUMINATION CONNECTOR

A502500	20 g
B502503	23 g
B502504	25 g



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

A502400	20 g
B502403	23 g
B502404	25 g



DORC EVA ILLUMINATION CONNECTOR

A502600	20 g
B502603	23 g
B502604	25 g



LIGHTMED LASER CONNECTOR



Protected by US Pat. 9,370,447, EP 2765943 and corresponding national patents.

ALCON STANDARD ILLUMINATION CONNECTOR

A50B200	20 g
B50B203	23 g
B50B204	25 g



19 GAUGE ILLUMINATION CONNECTOR

A50B100	20 g
B50B103	23 g
B50B104	25 g



ALCON OCUTOME ILLUMINATION CONNECTOR

A50B300	20 g
B50B303	23 g
B50B304	25 g



B&L STELLARIS PC ILLUMINATION CONNECTOR

A50B500	20 g
B50B503*	23 g
B50B504	25 g



B&L / STORZ MILLENIUM ILLUMINATION CONNECTOR

A50B400	20 g
B50B403	23 g
B50B404	25 g



DORC EVA ILLUMINATION CONNECTOR

A50B600	20 g
B50B603*	23 g
B50B604	25 g



* Not CE marked yet

CE
0413

PRODUCT INDEX

CODE	DESCRIPTION	PAGE	CODE	DESCRIPTION	PAGE	CODE	DESCRIPTION	PAGE
A202000	Bending 20g, Nidek thr.	3	A373000*	Laser 25g Angled, DORC	1	A523500	Illum. 20g Straight, DORC / Stellaris	7
A203000	Bending 20g, DORC	3	A374000	Laser 25g Angled, HGM	2	A523600	Illum. 20g Straight, DORC / DORC EVA	7
A207000	Bending 20g, Alcon/Zeiss	3	A375000	Laser 25g Angled, Nidek push	2	A524100	Illum. 20g Straight, HGM / 19 ga	12
A208000	Bending 20g, Iridex	3	A377000	Laser 25g Angled, Alcon/Zeiss	1	A524200	Illum. 20g Straight, HGM / Alcon Std.	12
A209000	Bending 20g, Quantel	3	A378000	Laser 25g Angled, Iridex	1	A524300	Illum. 20g Straight, HGM / Ocutome	12
A20A000	Bending 20g, Lumenis	3	A379000	Laser 25g Angled, Quantel	1	A524400	Illum. 20g Straight, HGM / Millenium	12
A20B000	Bending 20g, Lightmed	3	A37A000	Laser 25g Angled, Lumenis	1	A524500	Illum. 20g Straight, HGM / Stellaris	12
A212000	Laser 20g Straight, Nidek thr.	1	A37B000	Laser 25g Angled, Lightmed	2	A524600	Illum. 20g Straight, HGM / DORC EVA	12
A213000	Laser 20g Straight, DORC	1	A502100	Illum. bend. 20g, Nidek thr. / 19 ga	16	A525100	Illum. 20g Straight, Nidek push / 19 ga	12
A214000	Laser 20g Straight, HGM	2	A502200	Illum. bend. 20g, Nidek thr. / Alcon Std.	16	A525200	Illum. 20g Straight, Nidek push / Alcon Std.	12
A215000	Laser 20g Straight, Nidek push	2	A502300	Illum. bend. 20g, Nidek thr. / Ocutome	16	A525300	Illum. 20g Straight, Nidek push / Ocutome	12
A217000	Laser 20g Straight, Alcon/Zeiss	1	A502400	Illum. bend. 20g, Nidek thr. / Millenium	16	A525400	Illum. 20g Straight, Nidek push / Millenium	12
A218000	Laser 20g Straight, Iridex	1	A502500	Illum. bend. 20g, Nidek thr. / Stellaris	16	A525500	Illum. 20g Straight, Nidek push / Stellaris	12
A219000	Laser 20g Straight, Quantel	1	A502600	Illum. bend. 20g, Nidek thr. / DORC EVA	16	A525600	Illum. 20g Straight, Nidek push / DORC EVA	12
A21A000	Laser 20g Straight, Lumenis	1	A503100	Illum. bend. 20g, DORC / 19 ga	15	A527100	Illum. 20g Straight, Alcon/Zeiss / 19 ga	5
A21B000	Laser 20g Straight, Lightmed	2	A503200	Illum. bend. 20g, DORC / Alcon Std.	15	A527200	Illum. 20g Straight, Alcon/Zeiss / Alcon Std.	5
A232000	Laser 20g Angled, Nidek thr.	1	A503300	Illum. bend. 20g, DORC / Ocutome	15	A527300	Illum. 20g Straight, Alcon/Zeiss / Ocutome	5
A233000	Laser 20g Angled, DORC	1	A503400	Illum. bend. 20g, DORC / Millenium	15	A527400	Illum. 20g Straight, Alcon/Zeiss / Millenium	5
A234000	Laser 20g Angled, HGM	2	A503500	Illum. bend. 20g, DORC / Stellaris	15	A527500	Illum. 20g Straight, Alcon/Zeiss / Stellaris	5
A235000	Laser 20g Angled, Nidek push	2	A503600	Illum. bend. 20g, DORC / DORC EVA	15	A527600	Illum. 20g Straight, Alcon/Zeiss / DORC EVA	5
A237000	Laser 20g Angled, Alcon/Zeiss	1	A507100	Illum. bend. 20g, Alcon/Zeiss / 19 ga	14	A528100	Illum. 20g Straight, Iridex / 19 ga	6
A238000	Laser 20g Angled, Iridex	1	A507200	Illum. bend. 20g, Alcon/Zeiss / Alcon Std.	14	A528200	Illum. 20g Straight, Iridex / Alcon Std.	6
A239000	Laser 20g Angled, Quantel	1	A507300	Illum. bend. 20g, Alcon/Zeiss / Ocutome	14	A528300	Illum. 20g Straight, Iridex / Ocutome	6
A23A000	Laser 20g Angled, Lumenis	1	A507400	Illum. bend. 20g, Alcon/Zeiss / Millenium	14	A528400	Illum. 20g Straight, Iridex / Millenium	6
A23B000	Laser 20g Angled, Lightmed	2	A507500	Illum. bend. 20g, Alcon/Zeiss / Stellaris	14	A528500	Illum. 20g Straight, Iridex / Stellaris	6
A252000	Laser 20g Curved, Nidek thr.	1	A507600	Illum. bend. 20g, Alcon/Zeiss / DORC EVA	14	A528600	Illum. 20g Straight, Iridex / DORC EVA	6
A253000	Laser 20g Curved, DORC	1	A508100	Illum. bend. 20g, Iridex / 19 ga	14	A529100	Illum. 20g Straight, Quantel / 19 ga	8
A254000	Laser 20g Curved, HGM	2	A508200	Illum. bend. 20g, Iridex / Alcon Std.	14	A529200	Illum. 20g Straight, Quantel / Alcon Std.	8
A255000	Laser 20g Curved, Nidek push	2	A508300	Illum. bend. 20g, Iridex / Ocutome	14	A529300	Illum. 20g Straight, Quantel / Ocutome	8
A257000	Laser 20g Curved, Alcon/Zeiss	1	A508400	Illum. bend. 20g, Iridex / Millenium	14	A529400	Illum. 20g Straight, Quantel / Millenium	8
A258000	Laser 20g Curved, Iridex	1	A508500	Illum. bend. 20g, Iridex / Stellaris	14	A529500	Illum. 20g Straight, Quantel / Stellaris	8
A259000	Laser 20g Curved, Quantel	1	A508600	Illum. bend. 20g, Iridex / DORC EVA	14	A529600	Illum. 20g Straight, Quantel / DORC EVA	8
A25A000	Laser 20g Curved, Lumenis	1	A509100	Illum. bend. 20g, Quantel / 19 ga	15	A52A100	Illum. 20g Straight, Lumenis / 19 ga	9
A25B000	Laser 20g Curved, Lightmed	2	A509200	Illum. bend. 20g, Quantel / Alcon Std.	15	A52A200	Illum. 20g Straight, Lumenis / Alcon Std.	9
A302000	Bending 25g, Nidek thr.	3	A509300	Illum. bend. 20g, Quantel / Ocutome	15	A52A300	Illum. 20g Straight, Lumenis / Ocutome	9
A303000	Bending 25g, DORC	3	A509400	Illum. bend. 20g, Quantel / Millenium	15	A52A400	Illum. 20g Straight, Lumenis / Millenium	9
A307000	Bending 25g, Alcon/Zeiss	3	A509500	Illum. bend. 20g, Quantel / Stellaris	15	A52A500	Illum. 20g Straight, Lumenis / Stellaris	9
A308000	Bending 25g, Iridex	3	A509600	Illum. bend. 20g, Quantel / DORC EVA	15	A52A600	Illum. 20g Straight, Lumenis / DORC EVA	9
A309000	Bending 25g, Quantel	3	A50A100	Illum. bend. 20g, Lumenis / 19 ga	15	A52B100	Illum. 20g Straight, Lightmed / 19 ga	11
A30A000	Bending 25g, Lumenis	3	A50A200	Illum. bend. 20g, Lumenis / Alcon Std.	15	A52B200	Illum. 20g Straight, Lightmed / Alcon Std.	11
A30B000	Bending 25g, Lightmed	3	A50A300	Illum. bend. 20g, Lumenis / Ocutome	15	A52B300	Illum. 20g Straight, Lightmed / Ocutome	11
A322000	Laser 25g Straight, Nidek thr.	1	A50A400	Illum. bend. 20g, Lumenis / Millenium	15	A52B400	Illum. 20g Straight, Lightmed / Millenium	11
A323000	Laser 25g Straight, DORC	1	A50A500	Illum. bend. 20g, Lumenis / Stellaris	15	A52B500	Illum. 20g Straight, Lightmed / Stellaris	11
A324000	Laser 25g Straight, HGM	2	A50A600	Illum. bend. 20g, Lumenis / DORC EVA	15	A52B600	Illum. 20g Straight, Lightmed / DORC EVA	11
A325000	Laser 25g Straight, Nidek push	2	A50B100	Illum. bend. 20g, Lightmed / 19 ga	16	A542100	Illum. 20g Angled, Nidek thr. / 19 ga	10
A327000	Laser 25g Straight, Alcon/Zeiss	1	A50B200	Illum. bend. 20g, Lightmed / Alcon Std.	16	A542200	Illum. 20g Angled, Nidek thr. / Alcon Std.	10
A328000	Laser 25g Straight, Iridex	1	A50B300	Illum. bend. 20g, Lightmed / Ocutome	16	A542300	Illum. 20g Angled, Nidek thr. / Ocutome	10
A329000	Laser 25g Straight, Quantel	1	A50B400	Illum. bend. 20g, Lightmed / Millenium	16	A542400	Illum. 20g Angled, Nidek thr. / Millenium	10
A32A000	Laser 25g Straight, Lumenis	1	A50B500	Illum. bend. 20g, Lightmed / Stellaris	16	A542500	Illum. 20g Angled, Nidek thr. / Stellaris	10
A32B000	Laser 25g Straight, Lightmed	2	A50B600	Illum. bend. 20g, Lightmed / DORC EVA	16	A542600	Illum. 20g Angled, Nidek thr. / DORC EVA	10
A362000	Laser 25g Curved, Nidek thr.	1	A522100	Illum. 20g Straight, Nidek thr. / 19 ga	10	A543100	Illum. 20g Angled, DORC / 19 ga	7
A363000	Laser 25g Curved, DORC	1	A522200	Illum. 20g Straight, Nidek thr. / Alcon Std.	10	A543200	Illum. 20g Angled, DORC / Alcon Std.	7
A364000	Laser 25g Curved, HGM	2	A522300	Illum. 20g Straight, Nidek thr. / Ocutome	10	A543300	Illum. 20g Angled, DORC / Ocutome	7
A365000	Laser 25g Curved, Nidek push	2	A522400	Illum. 20g Straight, Nidek thr. / Millenium	10	A543400	Illum. 20g Angled, DORC / Millenium	7
A367000	Laser 25g Curved, Alcon/Zeiss	1	A522500	Illum. 20g Straight, Nidek thr. / Stellaris	10	A543500	Illum. 20g Angled, DORC / Stellaris	7
A368000	Laser 25g Curved, Iridex	1	A522600	Illum. 20g Straight, Nidek thr. / DORC EVA	10	A543600	Illum. 20g Angled, DORC / DORC EVA	7
A369000	Laser 25g Curved, Quantel	1	A523100	Illum. 20g Straight, DORC / 19 ga	7	A544100	Illum. 20g Angled, HGM / 19 ga	12
A36A000	Laser 25g Curved, Lumenis	1	A523200	Illum. 20g Straight, DORC / Alcon Std.	7	A544200	Illum. 20g Angled, HGM / Alcon Std.	12
A36B000	Laser 25g Curved, Lightmed	2	A523300	Illum. 20g Straight, DORC / Ocutome	7	A544300	Illum. 20g Angled, HGM / Ocutome	12
A372000	Laser 25g Angled, Nidek thr.	1	A523400	Illum. 20g Straight, DORC / Millenium	7	A544400	Illum. 20g Angled, HGM / Millenium	12

* Not CE marked yet

PRODUCT INDEX

CODE	DESCRIPTION	PAGE	CODE	DESCRIPTION	PAGE	CODE	DESCRIPTION	PAGE
A544500	Illum. 20g Angled, HGM / Stellaris	12	A565500	Illum. 20g Curved, Nidek push / Stellaris	12	AA08000	Bending 27g, Iridex	3
A544600	Illum. 20g Angled, HGM / DORC EVA	12	A565600	Illum. 20g Curved, Nidek push / DORC EVA	12	AA09000	Bending 27g, Quantel	3
A545100	Illum. 20g Angled, Nidek push / 19 ga	12	A567100	Illum. 20g Curved, Alcon/Zeiss / 19 ga	5	AA0A000	Bending 27g, Lumenis	3
A545200	Illum. 20g Angled, Nidek push / Alcon Std.	12	A567200	Illum. 20g Curved, Alcon/Zeiss / Alcon Std.	5	AA0B000	Bending 27g, Lightmed	3
A545300	Illum. 20g Angled, Nidek push / Ocutome	12	A567300	Illum. 20g Curved, Alcon/Zeiss / Ocutome	5	AA22000	Laser 27g Straight, Nidek thr.	1
A545400	Illum. 20g Angled, Nidek push / Millenium	12	A567400	Illum. 20g Curved, Alcon/Zeiss / Millenium	5	AA23000	Laser 27g Straight, DORC	1
A545500	Illum. 20g Angled, Nidek push / Stellaris	12	A567500	Illum. 20g Curved, Alcon/Zeiss / Stellaris	5	AA24000	Laser 27g Straight, HGM	2
A545600	Illum. 20g Angled, Nidek push / DORC EVA	12	A567600	Illum. 20g Curved, Alcon/Zeiss / DORC EVA	5	AA25000	Laser 27g Straight, Nidek push	2
A547100	Illum. 20g Angled, Alcon/Zeiss / 19 ga	5	A568100	Illum. 20g Curved, Iridex / 19 ga	6	AA27000	Laser 27g Straight, Alcon/Zeiss	1
A547200	Illum. 20g Angled, Alcon/Zeiss / Alcon Std.	5	A568200	Illum. 20g Curved, Iridex / Alcon Std.	6	AA28000	Laser 27g Straight, Iridex	1
A547300	Illum. 20g Angled, Alcon/Zeiss / Ocutome	5	A568300	Illum. 20g Curved, Iridex / Ocutome	6	AA29000	Laser 27g Straight, Quantel	1
A547400	Illum. 20g Angled, Alcon/Zeiss / Millenium	5	A568400	Illum. 20g Curved, Iridex / Millenium	6	AA2A000	Laser 27g Straight, Lumenis	1
A547500	Illum. 20g Angled, Alcon/Zeiss / Stellaris	5	A568500	Illum. 20g Curved, Iridex / Stellaris	6	AA2B000	Laser 27g Straight, Lightmed	2
A547600	Illum. 20g Angled, Alcon/Zeiss / DORC EVA	5	A568600	Illum. 20g Curved, Iridex / DORC EVA	6	AA62000	Laser 27g Curved, Nidek thr.	1
A548100	Illum. 20g Angled, Iridex / 19 ga	6	A569100	Illum. 20g Curved, Quantel / 19 ga	8	AA63000	Laser 27g Curved, DORC	1
A548200	Illum. 20g Angled, Iridex / Alcon Std.	6	A569200	Illum. 20g Curved, Quantel / Alcon Std.	8	AA64000	Laser 27g Curved, HGM	2
A548300	Illum. 20g Angled, Iridex / Ocutome	6	A569300	Illum. 20g Curved, Quantel / Ocutome	8	AA65000	Laser 27g Curved, Nidek push	2
A548400	Illum. 20g Angled, Iridex / Millenium	6	A569400	Illum. 20g Curved, Quantel / Millenium	8	AA67000	Laser 27g Curved, Alcon/Zeiss	1
A548500	Illum. 20g Angled, Iridex / Stellaris	6	A569500	Illum. 20g Curved, Quantel / Stellaris	8	AA68000	Laser 27g Curved, Iridex	1
A548600	Illum. 20g Angled, Iridex / DORC EVA	6	A569600	Illum. 20g Curved, Quantel / DORC EVA	8	AA69000	Laser 27g Curved, Quantel	1
A549100	Illum. 20g Angled, Quantel / 19 ga	8	A56A100	Illum. 20g Curved, Lumenis / 19 ga	9	AA6A000	Laser 27g Curved, Lumenis	1
A549200	Illum. 20g Angled, Quantel / Alcon Std.	8	A56A200	Illum. 20g Curved, Lumenis / Alcon Std.	9	AA6B000	Laser 27g Curved, Lightmed	2
A549300	Illum. 20g Angled, Quantel / Ocutome	8	A56A300	Illum. 20g Curved, Lumenis / Ocutome	9	B502103	Illum. bend. 23g, Nidek thr. / 19 ga	16
A549400	Illum. 20g Angled, Quantel / Millenium	8	A56A400	Illum. 20g Curved, Lumenis / Millenium	9	B502104	Illum. bend. 25g, Nidek thr. / 19 ga	16
A549500	Illum. 20g Angled, Quantel / Stellaris	8	A56A500	Illum. 20g Curved, Lumenis / Stellaris	9	B502203	Illum. bend. 23g, Nidek thr. / Alcon Std.	16
A549600	Illum. 20g Angled, Quantel / DORC EVA	8	A56A600	Illum. 20g Curved, Lumenis / DORC EVA	9	B502204	Illum. bend. 25g, Nidek thr. / Alcon Std.	16
A54A100	Illum. 20g Angled, Lumenis / 19 ga	9	A56B100	Illum. 20g Curved, Lightmed / 19 ga	11	B502303	Illum. bend. 23g, Nidek thr. / Ocutome	16
A54A200	Illum. 20g Angled, Lumenis / Alcon Std.	9	A56B200	Illum. 20g Curved, Lightmed / Alcon Std.	11	B502304	Illum. bend. 25g, Nidek thr. / Ocutome	16
A54A300	Illum. 20g Angled, Lumenis / Ocutome	9	A56B300	Illum. 20g Curved, Lightmed / Ocutome	11	B502403	Illum. bend. 23g, Nidek thr. / Millenium	16
A54A400	Illum. 20g Angled, Lumenis / Millenium	9	A56B400	Illum. 20g Curved, Lightmed / Millenium	11	B502404	Illum. bend. 25g, Nidek thr. / Millenium	16
A54A500	Illum. 20g Angled, Lumenis / Stellaris	9	A56B500	Illum. 20g Curved, Lightmed / Stellaris	11	B502503	Illum. bend. 23g, Nidek thr. / Stellaris	16
A54A600	Illum. 20g Angled, Lumenis / DORC EVA	9	A56B600	Illum. 20g Curved, Lightmed / DORC EVA	11	B502504	Illum. bend. 25g, Nidek thr. / Stellaris	16
A54B100	Illum. 20g Angled, Lightmed / 19 ga	11	A802000	Bending 23g, Nidek thr.	3	B502603	Illum. bend. 23g, Nidek thr. / DORC EVA	16
A54B200	Illum. 20g Angled, Lightmed / Alcon Std.	11	A803000	Bending 23g, DORC	3	B502604	Illum. bend. 25g, Nidek thr. / DORC EVA	16
A54B300	Illum. 20g Angled, Lightmed / Ocutome	11	A807000	Bending 23g, Alcon/Zeiss	3	B503103	Illum. bend. 23g, DORC / 19 ga	15
A54B400	Illum. 20g Angled, Lightmed / Millenium	11	A808000	Bending 23g, Iridex	3	B503104	Illum. bend. 25g, DORC / 19 ga	15
A54B500	Illum. 20g Angled, Lightmed / Stellaris	11	A809000	Bending 23g, Quantel	3	B503203	Illum. bend. 23g, DORC / Alcon Std.	15
A54B600	Illum. 20g Angled, Lightmed / DORC EVA	11	A80A000	Bending 23g, Lumenis	3	B503204	Illum. bend. 25g, DORC / Alcon Std.	15
A562100	Illum. 20g Curved, Nidek thr. / 19 ga	10	A80B000	Bending 23g, Lightmed	3	B503303	Illum. bend. 23g, DORC / Ocutome	15
A562200	Illum. 20g Curved, Nidek thr. / Alcon Std.	10	A822000	Laser 23g Straight, Nidek thr.	1	B503304	Illum. bend. 25g, DORC / Ocutome	15
A562300	Illum. 20g Curved, Nidek thr. / Ocutome	10	A823000	Laser 23g Straight, DORC	1	B503403	Illum. bend. 23g, DORC / Millenium	15
A562400	Illum. 20g Curved, Nidek thr. / Millenium	10	A824000	Laser 23g Straight, HGM	2	B503404	Illum. bend. 25g, DORC / Millenium	15
A562500	Illum. 20g Curved, Nidek thr. / Stellaris	10	A825000	Laser 23g Straight, Nidek push	2	B503503	Illum. bend. 23g, DORC / Stellaris	15
A562600	Illum. 20g Curved, Nidek thr. / DORC EVA	10	A827000	Laser 23g Straight, Alcon/Zeiss	1	B503504	Illum. bend. 25g, DORC / Stellaris	15
A563100	Illum. 20g Curved, DORC / 19 ga	7	A828000	Laser 23g Straight, Iridex	1	B503603	Illum. bend. 23g, DORC / DORC EVA	15
A563200	Illum. 20g Curved, DORC / Alcon Std.	7	A829000	Laser 23g Straight, Quantel	1	B503604	Illum. bend. 25g, DORC / DORC EVA	15
A563300	Illum. 20g Curved, DORC / Ocutome	7	A82A000	Laser 23g Straight, Lumenis	1	B507103	Illum. bend. 23g, Alcon/Zeiss / 19 ga	14
A563400	Illum. 20g Curved, DORC / Millenium	7	A82B000	Laser 23g Straight, Lightmed	2	B507104	Illum. bend. 25g, Alcon/Zeiss / 19 ga	14
A563500	Illum. 20g Curved, DORC / Stellaris	7	A872000	Laser 23g Angled, Nidek thr.	1	B507203	Illum. bend. 23g, Alcon/Zeiss / Alcon Std.	14
A563600	Illum. 20g Curved, DORC / DORC EVA	7	A873000	Laser 23g Angled, DORC	1	B507204	Illum. bend. 25g, Alcon/Zeiss / Alcon Std.	14
A564100	Illum. 20g Curved, HGM / 19 ga	12	A874000	Laser 23g Angled, HGM	2	B507205	Illum. bend. 27g, Alcon/Zeiss / Alcon Std.	14
A564200	Illum. 20g Curved, HGM / Alcon Std.	12	A875000	Laser 23g Angled, Nidek push	2	B507303	Illum. bend. 23g, Alcon/Zeiss / Ocutome	14
A564300	Illum. 20g Curved, HGM / Ocutome	12	A877000	Laser 23g Angled, Alcon/Zeiss	1	B507304	Illum. bend. 25g, Alcon/Zeiss / Ocutome	14
A564400	Illum. 20g Curved, HGM / Millenium	12	A878000	Laser 23g Angled, Iridex	1	B507305	Illum. bend. 27g, Alcon/Zeiss / Ocutome	14
A564500	Illum. 20g Curved, HGM / Stellaris	12	A879000	Laser 23g Angled, Quantel	1	B507403	Illum. bend. 23g, Alcon/Zeiss / Millenium	14
A564600	Illum. 20g Curved, HGM / DORC EVA	12	A87A000	Laser 23g Angled, Lumenis	1	B507404	Illum. bend. 25g, Alcon/Zeiss / Millenium	14
A565100	Illum. 20g Curved, Nidek push / 19 ga	12	A87B000	Laser 23g Angled, Lightmed	2	B507503	Illum. bend. 23g, Alcon/Zeiss / Stellaris	14
A565200	Illum. 20g Curved, Nidek push / Alcon Std.	12	AA02000	Bending 27g, Nidek thr.	3	B507504	Illum. bend. 25g, Alcon/Zeiss / Stellaris	14
A565300	Illum. 20g Curved, Nidek push / Ocutome	12	AA03000	Bending 27g, DORC	3	B507603	Illum. bend. 23g, Alcon/Zeiss / DORC EVA	14
A565400	Illum. 20g Curved, Nidek push / Millenium	12	AA07000	Bending 27g, Alcon/Zeiss	3	B507604	Illum. bend. 25g, Alcon/Zeiss / DORC EVA	14

PRODUCT INDEX

CODE	DESCRIPTION	PAGE	CODE	DESCRIPTION	PAGE	CODE	DESCRIPTION	PAGE
B508103	Illum. bend. 23g, Iridex / 19 ga	14	B522503	Illum. 23g Straight, Nidek thr. / Stellaris	10	B528303	Illum. 23g Straight, Iridex / Ocutoime	6
B508104	Illum. bend. 25g, Iridex / 19 ga	14	B522504	Illum. 25g Straight, Nidek thr. / Stellaris	10	B528304	Illum. 25g Straight, Iridex / Ocutoime	6
B508203	Illum. bend. 23g, Iridex / Alcon Std.	14	B522505	Illum. 27g Straight, Nidek thr. / Stellaris	10	B528305	Illum. 27g Straight, Iridex / Ocutoime	6
B508204	Illum. bend. 25g, Iridex / Alcon Std.	14	B522603	Illum. 23g Straight, Nidek thr. / DORC EVA	10	B528403	Illum. 23g Straight, Iridex / Millenium	6
B508303	Illum. bend. 23g, Iridex / Ocutoime	14	B522604	Illum. 25g Straight, Nidek thr. / DORC EVA	10	B528404	Illum. 25g Straight, Iridex / Millenium	6
B508304	Illum. bend. 25g, Iridex / Ocutoime	14	B522605	Illum. 27g Straight, Nidek thr. / DORC EVA	10	B528405	Illum. 27g Straight, Iridex / Millenium	6
B508403	Illum. bend. 23g, Iridex / Millenium	14	B523103	Illum. 23g Straight, DORC / 19 ga	7	B528503	Illum. 23g Straight, Iridex / Stellaris	6
B508404	Illum. bend. 25g, Iridex / Millenium	14	B523104	Illum. 25g Straight, DORC / 19 ga	7	B528504	Illum. 25g Straight, Iridex / Stellaris	6
B508503	Illum. bend. 23g, Iridex / Stellaris	14	B523105	Illum. 27g Straight, DORC / 19 ga	7	B528505	Illum. 27g Straight, Iridex / Stellaris	6
B508504	Illum. bend. 25g, Iridex / Stellaris	14	B523203	Illum. 23g Straight, DORC / Alcon Std.	7	B528603	Illum. 23g Straight, Iridex / DORC EVA	6
B508603	Illum. bend. 23g, Iridex / DORC EVA	14	B523204	Illum. 25g Straight, DORC / Alcon Std.	7	B528604	Illum. 25g Straight, Iridex / DORC EVA	6
B508604	Illum. bend. 25g, Iridex / DORC EVA	14	B523205	Illum. 27g Straight, DORC / Alcon Std.	7	B528605	Illum. 27g Straight, Iridex / DORC EVA	6
B509103	Illum. bend. 23g, Quantel / 19 ga	15	B523303	Illum. 23g Straight, DORC / Ocutoime	7	B529103	Illum. 23g Straight, Quantel / 19 ga	8
B509104	Illum. bend. 25g, Quantel / 19 ga	15	B523304	Illum. 25g Straight, DORC / Ocutoime	7	B529104	Illum. 25g Straight, Quantel / 19 ga	8
B509203	Illum. bend. 23g, Quantel / Alcon Std.	15	B523305	Illum. 27g Straight, DORC / Ocutoime	7	B529105	Illum. 27g Straight, Quantel / 19 ga	8
B509204	Illum. bend. 25g, Quantel / Alcon Std.	15	B523403	Illum. 23g Straight, DORC / Millenium	7	B529203	Illum. 23g Straight, Quantel / Alcon Std.	8
B509303	Illum. bend. 23g, Quantel / Ocutoime	15	B523404	Illum. 25g Straight, DORC / Millenium	7	B529204	Illum. 25g Straight, Quantel / Alcon Std.	8
B509304	Illum. bend. 25g, Quantel / Ocutoime	15	B523405	Illum. 27g Straight, DORC / Millenium	7	B529205	Illum. 27g Straight, Quantel / Alcon Std.	8
B509403	Illum. bend. 23g, Quantel / Millenium	15	B523503	Illum. 23g Straight, DORC / Stellaris	7	B529303	Illum. 23g Straight, Quantel / Ocutoime	8
B509404	Illum. bend. 25g, Quantel / Millenium	15	B523504	Illum. 25g Straight, DORC / Stellaris	7	B529304	Illum. 25g Straight, Quantel / Ocutoime	8
B509503	Illum. bend. 23g, Quantel / Stellaris	15	B523505	Illum. 27g Straight, DORC / Stellaris	7	B529305	Illum. 27g Straight, Quantel / Ocutoime	8
B509504	Illum. bend. 25g, Quantel / Stellaris	15	B523603	Illum. 23g Straight, DORC / DORC EVA	7	B529403	Illum. 23g Straight, Quantel / Millenium	8
B509603	Illum. bend. 23g, Quantel / DORC EVA	15	B523604	Illum. 25g Straight, DORC / DORC EVA	7	B529404	Illum. 25g Straight, Quantel / Millenium	8
B509604	Illum. bend. 25g, Quantel / DORC EVA	15	B523605	Illum. 27g Straight, DORC / DORC EVA	7	B529405	Illum. 27g Straight, Quantel / Millenium	8
B50A103	Illum. bend. 23g, Lumenis / 19 ga	15	B524103	Illum. 23g Straight, HGM / 19 ga	12	B529503	Illum. 23g Straight, Quantel / Stellaris	8
B50A104	Illum. bend. 25g, Lumenis / 19 ga	15	B524203	Illum. 23g Straight, HGM / Alcon Std.	12	B529504	Illum. 25g Straight, Quantel / Stellaris	8
B50A203	Illum. bend. 23g, Lumenis / Alcon Std.	15	B524303	Illum. 23g Straight, HGM / Ocutoime	12	B529505	Illum. 27g Straight, Quantel / Stellaris	8
B50A204	Illum. bend. 25g, Lumenis / Alcon Std.	15	B524403	Illum. 23g Straight, HGM / Millenium	12	B529603	Illum. 23g Straight, Quantel / DORC EVA	8
B50A303	Illum. bend. 23g, Lumenis / Ocutoime	15	B524503	Illum. 23g Straight, HGM / Stellaris	12	B529604	Illum. 25g Straight, Quantel / DORC EVA	8
B50A304	Illum. bend. 25g, Lumenis / Ocutoime	15	B524603	Illum. 23g Straight, HGM / DORC EVA	12	B529605	Illum. 27g Straight, Quantel / DORC EVA	8
B50A403	Illum. bend. 23g, Lumenis / Millenium	15	B525103	Illum. 23g Straight, Nidek push / 19 ga	12	B52A103	Illum. 23g Straight, Lumenis / 19 ga	9
B50A404	Illum. bend. 25g, Lumenis / Millenium	15	B525203	Illum. 23g Straight, Nidek push / Alcon Std.	12	B52A104	Illum. 25g Straight, Lumenis / 19 ga	9
B50A503	Illum. bend. 23g, Lumenis / Stellaris	15	B525303	Illum. 23g Straight, Nidek push / Ocutoime	12	B52A105	Illum. 27g Straight, Lumenis / 19 ga	9
B50A504	Illum. bend. 25g, Lumenis / Stellaris	15	B525403	Illum. 23g Straight, Nidek push / Millenium	12	B52A203	Illum. 23g Straight, Lumenis / Alcon Std.	9
B50A603	Illum. bend. 23g, Lumenis / DORC EVA	15	B525503	Illum. 23g Straight, Nidek push / Stellaris	12	B52A204	Illum. 25g Straight, Lumenis / Alcon Std.	9
B50A604	Illum. bend. 25g, Lumenis / DORC EVA	15	B525603	Illum. 23g Straight, Nidek push / DORC EVA	12	B52A205	Illum. 27g Straight, Lumenis / Alcon Std.	9
B50B103	Illum. bend. 23g, Lightmed / 19 ga	16	B527103	Illum. 23g Straight, Alcon/Zeiss / 19 ga	5	B52A303	Illum. 23g Straight, Lumenis / Ocutoime	9
B50B104	Illum. bend. 25g, Lightmed / 19 ga	16	B527104	Illum. 25g Straight, Alcon/Zeiss / 19 ga	5	B52A304	Illum. 25g Straight, Lumenis / Ocutoime	9
B50B203	Illum. bend. 23g, Lightmed / Alcon Std.	16	B527105	Illum. 27g Straight, Alcon/Zeiss / 19 ga	5	B52A305	Illum. 27g Straight, Lumenis / Ocutoime	9
B50B204	Illum. bend. 25g, Lightmed / Alcon Std.	16	B527203	Illum. 23g Straight, Alcon/Zeiss / Alcon Std.	5	B52A403	Illum. 23g Straight, Lumenis / Millenium	9
B50B303	Illum. bend. 23g, Lightmed / Ocutoime	16	B527204	Illum. 25g Straight, Alcon/Zeiss / Alcon Std.	5	B52A404	Illum. 25g Straight, Lumenis / Millenium	9
B50B304	Illum. bend. 25g, Lightmed / Ocutoime	16	B527205	Illum. 27g Straight, Alcon/Zeiss / Alcon Std.	5	B52A405	Illum. 27g Straight, Lumenis / Millenium	9
B50B403	Illum. bend. 23g, Lightmed / Millenium	16	B527303	Illum. 23g Straight, Alcon/Zeiss / Ocutoime	5	B52A503	Illum. 23g Straight, Lumenis / Stellaris	9
B50B404	Illum. bend. 25g, Lightmed / Millenium	16	B527304	Illum. 25g Straight, Alcon/Zeiss / Ocutoime	5	B52A504	Illum. 25g Straight, Lumenis / Stellaris	9
B50B503*	Illum. bend. 23g, Lightmed / Stellaris	16	B527305	Illum. 27g Straight, Alcon/Zeiss / Ocutoime	5	B52A505	Illum. 27g Straight, Lumenis / Stellaris	9
B50B504	Illum. bend. 25g, Lightmed / Stellaris	16	B527403	Illum. 23g Straight, Alcon/Zeiss / Millenium	5	B52A603	Illum. 23g Straight, Lumenis / DORC EVA	9
B50B603*	Illum. bend. 23g, Lightmed / DORC EVA	16	B527404	Illum. 25g Straight, Alcon/Zeiss / Millenium	5	B52A604	Illum. 25g Straight, Lumenis / DORC EVA	9
B50B604	Illum. bend. 25g, Lightmed / DORC EVA	16	B527405	Illum. 27g Straight, Alcon/Zeiss / Millenium	5	B52A605	Illum. 27g Straight, Lumenis / DORC EVA	9
B522103	Illum. 23g Straight, Nidek thr. / 19 ga	10	B527503	Illum. 23g Straight, Alcon/Zeiss / Stellaris	5	B52B103	Illum. 23g Straight, Lightmed / 19 ga	11
B522104	Illum. 25g Straight, Nidek thr. / 19 ga	10	B527504	Illum. 25g Straight, Alcon/Zeiss / Stellaris	5	B52B104	Illum. 25g Straight, Lightmed / 19 ga	11
B522105	Illum. 27g Straight, Nidek thr. / 19 ga	10	B527505	Illum. 27g Straight, Alcon/Zeiss / Stellaris	5	B52B105	Illum. 27g Straight, Lightmed / 19 ga	11
B522203	Illum. 23g Straight, Nidek thr. / Alcon Std.	10	B527603	Illum. 23g Straight, Alcon/Zeiss / DORC EVA	5	B52B203	Illum. 23g Straight, Lightmed / Alcon Std.	11
B522204	Illum. 25g Straight, Nidek thr. / Alcon Std.	10	B527604	Illum. 25g Straight, Alcon/Zeiss / DORC EVA	5	B52B204	Illum. 25g Straight, Lightmed / Alcon Std.	11
B522205	Illum. 27g Straight, Nidek thr. / Alcon Std.	10	B527605	Illum. 27g Straight, Alcon/Zeiss / DORC EVA	5	B52B205	Illum. 27g Straight, Lightmed / Alcon Std.	11
B522303	Illum. 23g Straight, Nidek thr. / Ocutoime	10	B528103	Illum. 23g Straight, Iridex / 19 ga	6	B52B303	Illum. 23g Straight, Lightmed / Ocutoime	11
B522304	Illum. 25g Straight, Nidek thr. / Ocutoime	10	B528104	Illum. 25g Straight, Iridex / 19 ga	6	B52B304	Illum. 25g Straight, Lightmed / Ocutoime	11
B522305	Illum. 27g Straight, Nidek thr. / Ocutoime	10	B528105	Illum. 27g Straight, Iridex / 19 ga	6	B52B305	Illum. 27g Straight, Lightmed / Ocutoime	11
B522403	Illum. 23g Straight, Nidek thr. / Millenium	10	B528203	Illum. 23g Straight, Iridex / Alcon Std.	6	B52B403	Illum. 23g Straight, Lightmed / Millenium	11
B522404	Illum. 25g Straight, Nidek thr. / Millenium	10	B528204	Illum. 25g Straight, Iridex / Alcon Std.	6	B52B404	Illum. 25g Straight, Lightmed / Millenium	11
B522405	Illum. 27g Straight, Nidek thr. / Millenium	10	B528205	Illum. 27g Straight, Iridex / Alcon Std.	6	B52B405	Illum. 27g Straight, Lightmed / Millenium	11

* Not CE marked yet

PRODUCT INDEX

CODE	DESCRIPTION	PAGE	CODE	DESCRIPTION	PAGE	CODE	DESCRIPTION	PAGE
B52B503	Illum. 23g Straight, Lightmed / Stellaris	11	B569404	Illum. 25g Curved, Quantel / Millennium	8	B577403	Illum. 23g Angled, Alcon/Zeiss / Millennium	5
B52B504	Illum. 25g Straight, Lightmed / Stellaris	11	B569405	Illum. 27g Curved, Quantel / Millennium	8	B577404	Illum. 25g Angled, Alcon/Zeiss / Millennium	5
B52B505	Illum. 27g Straight, Lightmed / Stellaris	11	B569504	Illum. 25g Curved, Quantel / Stellaris	8	B577503	Illum. 23g Angled, Alcon/Zeiss / Stellaris	5
B52B603	Illum. 23g Straight, Lightmed / DORC EVA	11	B569505	Illum. 27g Curved, Quantel / Stellaris	8	B577504	Illum. 25g Angled, Alcon/Zeiss / Stellaris	5
B52B604	Illum. 25g Straight, Lightmed / DORC EVA	11	B569604	Illum. 25g Curved, Quantel / DORC EVA	8	B577603	Illum. 23g Angled, Alcon/Zeiss / DORC EVA	5
B52B605	Illum. 27g Straight, Lightmed / DORC EVA	11	B569605	Illum. 27g Curved, Quantel / DORC EVA	8	B577604	Illum. 25g Angled, Alcon/Zeiss / DORC EVA	5
B562104	Illum. 25g Curved, Nidek thr. / 19 ga	10	B56A104	Illum. 25g Curved, Lumenis / 19 ga	9	B578103	Illum. 23g Angled, Iridex / 19 ga	6
B562105	Illum. 27g Curved, Nidek thr. / 19 ga	10	B56A105	Illum. 27g Curved, Lumenis / 19 ga	9	B578104	Illum. 25g Angled, Iridex / 19 ga	6
B562204	Illum. 25g Curved, Nidek thr. / Alcon Std.	10	B56A204	Illum. 25g Curved, Lumenis / Alcon Std.	9	B578203	Illum. 23g Angled, Iridex / Alcon Std.	6
B562205	Illum. 27g Curved, Nidek thr. / Alcon Std.	10	B56A205	Illum. 27g Curved, Lumenis / Alcon Std.	9	B578204	Illum. 25g Angled, Iridex / Alcon Std.	6
B562304	Illum. 25g Curved, Nidek thr. / Ocutoome	10	B56A304	Illum. 25g Curved, Lumenis / Ocutoome	9	B578303	Illum. 23g Angled, Iridex / Ocutoome	6
B562305	Illum. 27g Curved, Nidek thr. / Ocutoome	10	B56A305	Illum. 27g Curved, Lumenis / Ocutoome	9	B578304	Illum. 25g Angled, Iridex / Ocutoome	6
B562404	Illum. 25g Curved, Nidek thr. / Millennium	10	B56A404	Illum. 25g Curved, Lumenis / Millennium	9	B578403	Illum. 23g Angled, Iridex / Millennium	6
B562405	Illum. 27g Curved, Nidek thr. / Millennium	10	B56A405	Illum. 27g Curved, Lumenis / Millennium	9	B578404	Illum. 25g Angled, Iridex / Millennium	6
B562504	Illum. 25g Curved, Nidek thr. / Stellaris	10	B56A504	Illum. 25g Curved, Lumenis / Stellaris	9	B578503	Illum. 23g Angled, Iridex / Stellaris	6
B562505	Illum. 27g Curved, Nidek thr. / Stellaris	10	B56A505	Illum. 27g Curved, Lumenis / Stellaris	9	B578504	Illum. 25g Angled, Iridex / Stellaris	6
B562604	Illum. 25g Curved, Nidek thr. / DORC EVA	10	B56A604	Illum. 25g Curved, Lumenis / DORC EVA	9	B578603	Illum. 23g Angled, Iridex / DORC EVA	6
B562605	Illum. 27g Curved, Nidek thr. / DORC EVA	10	B56A605	Illum. 27g Curved, Lumenis / DORC EVA	9	B578604	Illum. 25g Angled, Iridex / DORC EVA	6
B563104	Illum. 25g Curved, DORC / 19 ga	7	B56B104	Illum. 25g Curved, Lightmed / 19 ga	11	B579103	Illum. 23g Angled, Quantel / 19 ga	8
B563105	Illum. 27g Curved, DORC / 19 ga	7	B56B105	Illum. 27g Curved, Lightmed / 19 ga	11	B579104	Illum. 25g Angled, Quantel / 19 ga	8
B563204	Illum. 25g Curved, DORC / Alcon Std.	7	B56B204	Illum. 25g Curved, Lightmed / Alcon Std.	11	B579203	Illum. 23g Angled, Quantel / Alcon Std.	8
B563205	Illum. 27g Curved, DORC / Alcon Std.	7	B56B205	Illum. 27g Curved, Lightmed / Alcon Std.	11	B579204	Illum. 25g Angled, Quantel / Alcon Std.	8
B563304	Illum. 25g Curved, DORC / Ocutoome	7	B56B304	Illum. 25g Curved, Lightmed / Ocutoome	11	B579303	Illum. 23g Angled, Quantel / Ocutoome	8
B563305	Illum. 27g Curved, DORC / Ocutoome	7	B56B305	Illum. 27g Curved, Lightmed / Ocutoome	11	B579304	Illum. 25g Angled, Quantel / Ocutoome	8
B563404	Illum. 25g Curved, DORC / Millennium	7	B56B404	Illum. 25g Curved, Lightmed / Millennium	11	B579403	Illum. 23g Angled, Quantel / Millennium	8
B563405	Illum. 27g Curved, DORC / Millennium	7	B56B405	Illum. 27g Curved, Lightmed / Millennium	11	B579404	Illum. 25g Angled, Quantel / Millennium	8
B563504	Illum. 25g Curved, DORC / Stellaris	7	B56B504	Illum. 25g Curved, Lightmed / Stellaris	11	B579503	Illum. 23g Angled, Quantel / Stellaris	8
B563505	Illum. 27g Curved, DORC / Stellaris	7	B56B505	Illum. 27g Curved, Lightmed / Stellaris	11	B579504	Illum. 25g Angled, Quantel / Stellaris	8
B563604	Illum. 25g Curved, DORC / DORC EVA	7	B56B604	Illum. 25g Curved, Lightmed / DORC EVA	11	B579603	Illum. 23g Angled, Quantel / DORC EVA	8
B563605	Illum. 27g Curved, DORC / DORC EVA	7	B56B605	Illum. 27g Curved, Lightmed / DORC EVA	11	B579604	Illum. 25g Angled, Quantel / DORC EVA	8
B567104	Illum. 25g Curved, Alcon/Zeiss / 19 ga	5	B572103	Illum. 23g Angled, Nidek thr. / 19 ga	10	B57A103	Illum. 23g Angled, Lumenis / 19 ga	9
B567105	Illum. 27g Curved, Alcon/Zeiss / 19 ga	5	B572104	Illum. 25g Angled, Nidek thr. / 19 ga	10	B57A104	Illum. 25g Angled, Lumenis / 19 ga	9
B567204	Illum. 25g Curved, Alcon/Zeiss / Alcon Std.	5	B572203	Illum. 23g Angled, Nidek thr. / Alcon Std.	10	B57A203	Illum. 23g Angled, Lumenis / Alcon Std.	9
B567205	Illum. 27g Curved, Alcon/Zeiss / Alcon Std.	5	B572204	Illum. 25g Angled, Nidek thr. / Alcon Std.	10	B57A204	Illum. 25g Angled, Lumenis / Alcon Std.	9
B567304	Illum. 25g Curved, Alcon/Zeiss / Ocutoome	5	B572303	Illum. 23g Angled, Nidek thr. / Ocutoome	10	B57A303*	Illum. 23g Angled, Lumenis / Ocutoome	9
B567305	Illum. 27g Curved, Alcon/Zeiss / Ocutoome	5	B572304	Illum. 25g Angled, Nidek thr. / Ocutoome	10	B57A304	Illum. 25g Angled, Lumenis / Ocutoome	9
B567404	Illum. 25g Curved, Alcon/Zeiss / Millennium	5	B572403	Illum. 23g Angled, Nidek thr. / Millennium	10	B57A403	Illum. 23g Angled, Lumenis / Millennium	9
B567405	Illum. 27g Curved, Alcon/Zeiss / Millennium	5	B572404	Illum. 25g Angled, Nidek thr. / Millennium	10	B57A404	Illum. 25g Angled, Lumenis / Millennium	9
B567504	Illum. 25g Curved, Alcon/Zeiss / Stellaris	5	B572503	Illum. 23g Angled, Nidek thr. / Stellaris	10	B57A503	Illum. 23g Angled, Lumenis / Stellaris	9
B567505	Illum. 27g Curved, Alcon/Zeiss / Stellaris	5	B572504	Illum. 25g Angled, Nidek thr. / Stellaris	10	B57A504	Illum. 25g Angled, Lumenis / Stellaris	9
B567604	Illum. 25g Curved, Alcon/Zeiss / DORC EVA	5	B572603	Illum. 23g Angled, Nidek thr. / DORC EVA	10	B57A603	Illum. 23g Angled, Lumenis / DORC EVA	9
B567605	Illum. 27g Curved, Alcon/Zeiss / DORC EVA	5	B572604	Illum. 25g Angled, Nidek thr. / DORC EVA	10	B57A604	Illum. 25g Angled, Lumenis / DORC EVA	9
B568104	Illum. 25g Curved, Iridex / 19 ga	6	B573103	Illum. 23g Angled, DORC / 19 ga	7	B57B103	Illum. 23g Angled, Lightmed / 19 ga	11
B568105	Illum. 27g Curved, Iridex / 19 ga	6	B573104	Illum. 25g Angled, DORC / 19 ga	7	B57B104	Illum. 25g Angled, Lightmed / 19 ga	11
B568204	Illum. 25g Curved, Iridex / Alcon Std.	6	B573203	Illum. 23g Angled, DORC / Alcon Std.	7	B57B203	Illum. 23g Angled, Lightmed / Alcon Std.	11
B568205	Illum. 27g Curved, Iridex / Alcon Std.	6	B573204	Illum. 25g Angled, DORC / Alcon Std.	7	B57B204	Illum. 25g Angled, Lightmed / Alcon Std.	11
B568304	Illum. 25g Curved, Iridex / Ocutoome	6	B573303	Illum. 23g Angled, DORC / Ocutoome	7	B57B303*	Illum. 23g Angled, Lightmed / Ocutoome	11
B568305	Illum. 27g Curved, Iridex / Ocutoome	6	B573304	Illum. 25g Angled, DORC / Ocutoome	7	B57B304	Illum. 25g Angled, Lightmed / Ocutoome	11
B568404	Illum. 25g Curved, Iridex / Millennium	6	B573403	Illum. 23g Angled, DORC / Millennium	7	B57B403	Illum. 23g Angled, Lightmed / Millennium	11
B568405	Illum. 27g Curved, Iridex / Millennium	6	B573404	Illum. 25g Angled, DORC / Millennium	7	B57B404*	Illum. 25g Angled, Lightmed / Millennium	11
B568504	Illum. 25g Curved, Iridex / Stellaris	6	B573503	Illum. 23g Angled, DORC / Stellaris	7	B57B503	Illum. 23g Angled, Lightmed / Stellaris	11
B568505	Illum. 27g Curved, Iridex / Stellaris	6	B573504	Illum. 25g Angled, DORC / Stellaris	7	B57B504	Illum. 25g Angled, Lightmed / Stellaris	11
B568604	Illum. 25g Curved, Iridex / DORC EVA	6	B573603	Illum. 23g Angled, DORC / DORC EVA	7	B57B603	Illum. 23g Angled, Lightmed / DORC EVA	11
B568605	Illum. 27g Curved, Iridex / DORC EVA	6	B573604	Illum. 25g Angled, DORC / DORC EVA	7	B57B604	Illum. 25g Angled, Lightmed / DORC EVA	11
B569104	Illum. 25g Curved, Quantel / 19 ga	8	B577103	Illum. 23g Angled, Alcon/Zeiss / 19 ga	5			
B569105	Illum. 27g Curved, Quantel / 19 ga	8	B577104	Illum. 25g Angled, Alcon/Zeiss / 19 ga	5			
B569204	Illum. 25g Curved, Quantel / Alcon Std.	8	B577203	Illum. 23g Angled, Alcon/Zeiss / Alcon Std.	5			
B569205	Illum. 27g Curved, Quantel / Alcon Std.	8	B577204	Illum. 25g Angled, Alcon/Zeiss / Alcon Std.	5			
B569304	Illum. 25g Curved, Quantel / Ocutoome	8	B577303	Illum. 23g Angled, Alcon/Zeiss / Ocutoome	5			
B569305	Illum. 27g Curved, Quantel / Ocutoome	8	B577304	Illum. 25g Angled, Alcon/Zeiss / Ocutoome	5			

* Not CE marked yet



Aspiration add-on functionality is available upon request.

A flexible plastic jacket covers and protects the 8ft long fiber.

Homogeneous energy distribution as well as round and consistently small laser spot allow the tip to stay a few millimeters away from the retina.

Ergonomic and non-slip handpiece design.



Different connectors that suit most laser and illumination machines.

Sterilized by Ethylene Oxide, double pouched, single-use, latex-free.

Available in Dia 0.41 mm (27ga), Dia 0.51 mm (25ga), Dia 0.64 mm (23ga), and Dia 0.91 mm (20g). Length 2.4 m (8 ft). Box of 10.

Do not re-sterilize. Do not reuse.

Material fatigue decreases performance.



WARNING

Follow the recommendations of the laser and illumination source manufacturers (e.g.: calibration of the laser source as well as the need for protective eye wear for the corresponding wavelength of the laser model used).

WARNING IN CASE OF ILLUMINATION

Due to the potentially hazardous effects of phototoxicity to the retina from prolonged endoillumination, this device should be used with an illumination source with filters that eliminate UV radiation (<400nm) and whenever possible, with filters that eliminate short wavelength blue light (<420nm). Care should be taken not to focus this device on a single point of the retina for prolonged, unnecessary periods of time.

CAUTION: FEDERAL (U.S.A.) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN

CONTACT INFORMATION:

+1 770 777 6613

+1 877 4 25 GAUGE (4284)

+1 678 623 3765 (FAX)

info@ophthalmed.com

1050 Northfield Court, Suite 280
Roswell, GA 30076 USA

DISTRIBUTED BY:



OPHTHALMEDTM

MADE IN THE USA